

PROGRAM DESCRIPTION GUIDE

A. IDENTIFICATION

Program Name - Grodins (Respiratory Control Model)  
Programmer's Name - Marks, Archer  
Programmer Contact - V. J. Marks or G. T. Archer, GE/TSS, Houston  
Date of Issue - 9/6/74

B. GENERAL DESCRIPTION

The purpose of this model is to illustrate the transient and steady-state responses of the respiratory control system, for variations in volumetric fractions of inspired gases and special system parameters. The program contains the capability to change workload.

The program is based on Grodins' respiratory control model and can be envisioned as a feedback control system comprised of a "plant" (the controlled system) and the regulating component (controlling system). The controlled system is partitioned into 3 compartments corresponding to lungs, brain, and tissue with a fluid interconnecting patch representing the blood.

C. USAGE AND RESTRICTIONS

Machine, Operating System, and  
Compiler Required - Univac 1110, EXEC 8, Fortran  
Peripheral Equipment Required - Printer, Card Reader, Graphic Terminal  
Approximate Memory Required - 18671

D. PARTICULAR DESCRIPTION

See TIR 741-MED-3047

E. DESCRIPTION OF INPUT

A Univac 1110 file (GRØDIN) contains the source and relocatables of all the subroutines, the executable program, Batch and Demand data files (GRØDIN.GRØDAT and GRØDIN.GRØDATF), and Remote Batch run stream (GRØDIN.RUNB). Runs may require a data file and/or the Remote Batch run stream file to be modified. Since GRØDIN is not protected, the user should copy

GRØDIN. into another file, and make his modifications to this other file. Examples to follow should explain.

1. Batch -

The print output will be generated at the 1110 (onsite).

a. Onsite Batch

The run will be in card deck form, and initiated at the 1110. A JSC form 588A should be submitted with the card deck.

(1) Job Stream Cards (Start Col. 1)

```
@RUN,/R DBHØLD,7007-Q509-C,DB6-G03432,TT,PPP
  (TT=mins.run time, PPP=pages output, NAME=user
  name)
@ASG,AX GRØDIN.
@CØPY GRØDIN.,TFF$.
@FREE GRØDIN.
@XQT
  data cards
@FIN
```

Col.61  
NAME SKYLAB

(2) Data Cards (See Appendix A for example).

<u>Card No.</u>	<u>Column</u>	<u>Format</u>	<u>Description</u>
1	1-4	A4	Blank if Batch Mode. TTYb if Demand Mode.
2-49			Cards 2-49 are initializing cards for 48 variables (See Appendix A).
	6-20	F15.0	Variable value.
	26-33	2A4	Variable symbol (this is not used in the program).
50-N			Cards 50-N are workload cards.
	1-6	F6.2	Workload (watts).
	10-15	F6.2	*Run time (mins) for workload
N+1			*Workload card with time = 0.

\*Length of computer run is determined when variable on card 16 (maximum length of computer run) is reached, or a workload card with 0 run time is read.

b. Remote Batch

The run is initiated from a teletype. The run stream and the data are in Univac 1110 files GRØDIN.RUNB and GRØDIN.GRØDAT. The data file is the same format as the Onsite Batch data cards. The following example shows copying GRØDIN. into another file

FILNEW.), modifying the run stream and data files, copying the run stream into another file (FILRUN.), then executing the run via a teletype.

<u>TTY Input</u>	<u>Comments</u>
@ASG,CP FILNEW.,F2	Catalog a new file named FILNEW.
@ASG,CP FILRUN.,F2	Catalog a new file to put run stream. Run stream must be in separate file to execute.
@COPY GRØDIN.,TPF\$.	Copy GRØDIN. into temporary work file.
@ED,L RUNB,RUNB	To edit run stream.

The run stream is the following format:

@RUN,/R DBHOLD,7007-Q509-C,DB6-G03432,TT,PPP (TT= mins. run time, PPP = pages output) @ASG,AX FILNEW. @COPY FILNEW.,TPF\$. @FREE FILNEW. @XQT @ADD GRØDAT	↙ Col. 61 NAME NODECK
---	--------------------------

The user might want to change the run time or output pages in this run stream.

@ED,L GRØDAT,GRØDAT	To edit the data file.
The data file is in the same format as Onsite Batch data cards. The user can modify this base data file for his particular run.	
@COPY TPF\$.,FILNEW.	Copy temporary work file into FILNEW..
@COPY TPF\$., FILRUN.	To get run stream into FILRUN..
@FREE FILRUN.	
@START FILRUN.RUN	Start Remote Batch execution.

NOTE: After above TTY input, FILNEW. contains the latest run stream and data files. FILNEW. is not a permanent file, but it might be several days before NASA deletes it. The next time the user wants to modify these files, he may substitute the first three TTY inputs with @COPY FILNEW.,TPF\$. until FILNEW. is deleted.

## 2. Demand (Time-Share) -

The run is initiated from a teletype, and the printout will be to the teletype. A data file (GRØDIN.GRØDATT) is required to initialize the 49 input variables. Workloads are then input by the user via responding to questions asked by the program. Options are available so all workloads can be input before the program executes, or intermediate results may be requested before continued input. Note following example (Appendix B also contains an example):

### TTY Input

### Comments

@COPY GRODIN.,TPF\$.

Copy GRØDIN. into temporary work file.

@ED,L GRODATT,GRODATT

To edit Demand Data File.

The Demand Data File is in the same format as Batch Data File except the first card record contains TTY in first three columns, and the 49th card record ends the data file (no workload cards). Card record 16 (maximum length of computer run) and record 40 (time increment for printout) are ignored by the program in the Demand Mode. The user can modify this data file for his particular run.

@XQT

Execute program.

ADD DATA...

Program request.

@ADD GRØDATT

User answer.

The 48 input variables will then be printed out.

INPUT WORK CARDS...

WØRK = WØRK LØAD (WATTS)...

MINS = TIME FØR WØRK LØAD...

PRINT = TIME INCREMENT (MINS) FOR PRINTØUT

EXEC...

MØRE = INPUT MØRE BEFORE EXEC...

RUN = EXEC. WITH ABOVE, THEN CAN INPUT AGAIN...

STØP = EXEC. WITH ABOVE, THEN STØP...

BACK = ERASE PREVIOUS WØRK RECØRD...

Program outputs an explanation of parameters on work card records.

<u>WØRK</u>	<u>MINS</u>	<u>PRINT</u>	<u>EXEC</u>	(F6.2,IX,F6.2,IX,F6.2,IX,A4)
0.	.2	.1	MØRE	
.00	.20	.10	MØRE	

Prog.asks for work record  
User answer.  
Prog. outputs user input.

<u>WØRK</u>	<u>MINS</u>	<u>PRINT</u>	<u>EXEC</u>	(F6.2,IX,F6.2,IX,F6.2,IX,A4)
100.	10.	.2	RUN	
100.00	10.00	.20	RUN	

Prog.asks for work record  
User answer.  
Prog.outputs user input.

The program executes and outputs the results for the preceding input, then (since EXEC was RUN) the program will ask for additional work records.

WORK    MINS    PRINT    EXEC (F6.2,IX,F6.2,IX,F6.2,IX,A4)...Prog.asks for  
work record  
Etc.

Until the program finishes executing a work record with an EXEC  
parameter = STOP.

#### F. DESCRIPTION OF OUTPUT

(See Appendix B for Demand example)  
(See Appendix C for Batch example)  
(See Appendix D for normal values)

TIME (mins)

ALVEOLAR, ARTERIAL, BRAIN, TISSUE, V BRAIN, and V TISSUE  
volumetric fractions of CO<sub>2</sub>, O<sub>2</sub>, and N<sub>2</sub>.

ALVEOLAR, ARTERIAL, BRAIN, TISSUE, CSF, V BRAIN, and V TISSUE  
partial pressures of CO<sub>2</sub>, O<sub>2</sub>, and N<sub>2</sub>. (mm Hg)

ALVEOLAR, BRAIN, TISSUE, and CSF derivatives of the partial pressures  
of CO<sub>2</sub>, O<sub>2</sub>, and N<sub>2</sub>. (mm Hg)

ARTERIAL, BRAIN, TISSUE, CSF, V BRAIN, and V TISSUE  
hydrogen ion (H<sup>+</sup>) concentrations (nanomoles) and pH.

ARTERIAL, V BRAIN, and V TISSUE concentrations of HbO<sub>2</sub> (oxyhemoglobin),  
(liters O<sub>2</sub> - STPD)

ALVEOLAR RQ

RQ DIFF

TRANSPORT TIMES (mins)

AB = Lung to brain  
VB = Brain to lung  
VT = Tissue to lung  
AT = Lung to tissue  
AC = Lung to carotid body

VI = Inspired ventilation (liters/min)

VE = Expired ventilation (liters/min)

Q = Cardiac output (liters/min)

FB = Brain blood flow (liters/min)

DERIVATIVES of Q and FB

RESP FREQ (breaths/min)

MINUTE VOLUME (liters/min)

DEAD SPACE VENTILATION (liters/min)

HEART RATE (beats/min)

ARTERIAL VENOUS  $O_2$  DIFFERENCE (liters  $O_2$ /liters blood)

DEAD SPACE VOLUME (liters)

CHANGE IN WORK LOAD (watts)

METABOLIC RATE CHANGE OF  $CO_2$  AND  $O_2$  CONSUMPTION IN TISSUE

At termination of a run, the final values for the 2nd thru 15th input variables (see Appendix A) are printed out. If the 38th input variable (constant involved in controller equation)  $\leq .00001$ , it is recalculated and printed at the termination of the run.

#### G. INTERNAL CHECKS AND EXITS

##### 1. Batch

The run terminates when time becomes greater than the 16th input variable (maximum length of computer run), or when a workload card with 0 time is read.

##### 2. Demand

Up to 50 work records can be input before letting the program execute. If this limit is reached, the program will execute using an EXEC parameter = RUN. The EXEC parameter is checked for validity. The run terminates when a work record with an EXEC parameter = STOP is finished executing.

#### H. INDEPENDENT SUBROUTINES

See Appendix E for listing of all subroutines.

#### I. SYSTEM SUBROUTINES

No special system subroutines required. Library File TEK. is required to supply terminal graphic subroutines.

J. COMPLETION OR FINAL CHECKOUT DATE

9/6/74

Appendix A - Data Card Example  
Appendix B - Demand Example  
Appendix C - Batch Example  
Appendix D - Normal Values  
Appendix E - Program listing

APPENDIX A  
INITIALIZING CARDS FOR  
INPUT VARIABLES



<u>Card No.</u>	<u>Variable</u>	<u>Normal Initial Value</u>	<u>Symbol</u>	<u>Description</u>
1	TTY			Blank for Batch Mode. TTY in 1st 3 cols. if Demand Mode.
2	C(1)	.17827	FA(CO2)	Alveolar gas fractions (dry), volumetric fraction of gas, dimensionless
3	C(2)	.53459	FA(O2)	
4	C(3)	.28714	FA(N2)	
5	C(4)	.64121	CB(CO2)	Concentration of gas in brain, liters (STPD)/liter brain.
6	C(5)	.00116	CB(O2)	
7	C(6)	.00105	CB(N2)	
8	C(7)	.61553	CT(CO2)	Concentration of gas in tissue compartment. Liters (STPD)/ liter tissue
9	C(8)	.00147	CT(O2)	
10	C(9)	.00105	CT(N2)	
11	C(10)	6.00000	Q	Cardiac output blood flow, liters/min. Cerebral blood flow, liters/min.
12	C(11)	.74913	QB	
13	C(12)	48.17427	PCSF(CO2)	Partial pressure of gas in cerebrospinal fluid compart- ment, mm Hg.
14	C(13)	36.69498	PCSF(O2)	
15	C(14)	61.17176	PCSF(N2)	
16	C(15)	40.0000	TMAX	Maximum length of computer run, min. A workload card with time=0 will also end a Batch computer run. Demand Mode ignores this variable.
17	C(16)	0.0000	CENT SENS PT	Central sensitivity partition. Weighting of the H <sup>+</sup> concentra- tion in CSF with that of venous blood in the brain. With C(16) =0, zero weight is given to venous blood at level of the brain and a weight of one is given to H <sup>+</sup> concentration in CSF.
18	C(17)	.2000	HB	Blood oxygen capacity, liters (STPD)/liter blood Time constants for cardiac output response (R1) and cere- bral blood flow response (R2) for changes in blood chemical composition.
19	C(18)	.1000	R1	
20	C(19)	.1000	R2	
21	C(20)	1.1380	CNT SENS COF	Controller sensitivity weight- ings; i.e.,
22	C(21)	1.1540	CRTD BDV SCF	

Card No.	Variable	Normal Initial Value	Symbol
----------	----------	----------------------	--------

$$V_1 = 1.138C_{CSF(H^+)} + 1.1540 \\ (t - \tau_{ao}) + \text{TERM} - V_{1(N)}$$

where

$\tau_{ao}$  = Blood transport delay from  
lung to carotid body,

$V_1$  defined in C(37), and

TERM = function of  $F_{A(O_2)}$ .

23	C(22)	3.0000	KL	Volumes of lung (alveoli), brain,
24	C(23)	1.0000	KB	and tissue compartments, liters.
25	C(24)	39.0000	KT	
26	C(25)	.0500	MRB (CO2)	Metabolic rates by brain,
27	C(26)	.0500	MRB (O2)	liters (STPD)/min.
28	C(27)	81.9900	D (CO2)	Diffusion coefficient for gas
29	C(28)	4.3610	D (O2)	across "blood-brain", liters(10) <sup>-7</sup>
30	C(29)	2.5240	D (N2)	(STPD)/min. per mm Hg.
31	C(30)	260.0000	B	Barometric pressure, mm Hg.
32	C(31)	.0192	F1 (CO2)	Volumetric fraction of gas (dry
33	C(32)	.7000	F1 (O2)	inspired), dimensionless
34	C(33)	.2808	F1 (N2)	
35	C(34)	.1000	KCSF	Volume of cerebrospinal fluid,
				liters
36	C(35)	.0000	T	Initial time.
37	C(36)	.0078125	H	Size of computer time step, min.
38	C(37)	87.5500	V1(N)	Constant that is involved in the controller equation (See C(21)). Determines the normal level of alveolar ventilation so that $P_{A(CO_2)} \approx 40.0$ at rest, breathing air at sea level. When the con- troller sensitivity weightings are changed V1(N) should be altered accordingly.
39	C(38)	5.3900	V1 (SS)	Value used for normal resting alveolar ventilation. This is not used in the program if V1(N) is known.

<u>Card No.</u>	<u>Variable</u>	<u>Normal Initial Value</u>	<u>Symbol</u>	<u>Description</u>
40	C(39)	.2500	PRINT AL TIM	Output printed in these time increments (mins). Demand Mode ignores this variable.
41	C(40)	0.0000	UNKNOWN	Importance related to C (39), but doesn't seem to be of any real significance.
42	BC(1)	.5470	BHC03 Blood	Standard bicarbonate content, liters CO <sub>2</sub> (STPD)/liter X, 37° where
43	BC(2)	.5850		
44	BC(3)	.5850		
45	BC(4)	.5850		
				X = Blood, brain, tissue, CSF.
46	RMT(1)	.1820	RMT(CO2)	Metabolic rates by tissue, liters (STPD)/min.
47	RMT(2)	.2150	RMT(O2)	
48	DJ(1)	.0000	DJ1	Used in performing Dejours experiment (not utilized in present runs). Brief description of Dejours work relating O <sub>2</sub> and CO <sub>2</sub> threshold effects is given in Grodins' paper.
49	DJ(2)	.0000	DJ2	

In Demand Mode, the 49th card record ends the data file.

If for Batch Mode, cards 50 thru N are workload cards.

Example:

<u>Card No.</u>	<u>Variable</u>	<u>Columns</u>	<u>Format</u>	<u>Description</u>
50-N	WORK2	1-6	F6.2	Workload (watts).
	DURAT	10-15	F6.2	Run time for workload (mins).

The last workload card should have Run Time = 0. Length of computer run (Batch Mode) will be when C(15) is reached, or a workload card with 0 run time is read.

## APPENDIX B

### EXAMPLE OF DEMAND INPUT/OUTPUT

@XQT

# GRODINS: RESPIRATORY CONTROL MODEL

ADD DATA...  
 >@ADD GRODATT

## \*RESPIRATORY CHEMOSTAT -- INPUT DATA\*

1	.1783	.5346	.2871	.6412	.0012
6	.0011	.6155	.0015	.0011	6.0000
11	.7491	48.1742	36.6949	61.1717	40.0000
16	.0000	.2000	.1000	.1000	1.1380
21	1.1540	3.0000	1.0000	39.0000	.0500
26	.0500	81.9900	4.3610	2.5240	260.0000
31	.0192	.7000	.2800	.1000	.0000
36	.0078	87.5500	5.3900	.2500	.0000
41	.5470	.5850	.5850	.5850	
45	.1820	.2150	.0000	.0000	

DO YOU WANT GRAPHIC INSTEAD OF TABULAR OUTPUT? (Y/N)  
 N

INPUT WORK CARDS...

WORK= WORK LOAD(WATTS)...

MINS= TIME FOR WORK LOAD...

PRINT= TIME INCRIMENT(MINS)FOR PRINTOUT...

EXEC

MORE= INPUT MORE BEFORE EXEC...

RUN = EXEC WITH ABOVE THEN CAN INPUT AGAIN...

STOP= EXEC WITH ABOVE THEN STOP...

BACK= ERASE PREVIOUS WORK RECORD...

WORK MINS PRINT EXEC (F6.2,1X,F6.2,1X,F6.2,1X,A4)...

0. 1. .5 STOP  
 .00 1.00 .50STOP

\*\*\*\*\*

WORK LOAD CHG.( .00WATTS FOR 1.00MINS) AT .0000MINS

TIME	.0000MINS	ALV RQ	.8782	RQ DIFF	-.0027			
	ALVEOLAR	ARTERIAL	BRAIN	TISSUE	CSF	U	BRAIN	U
							TISSUE	
CO2	.1783	.5652	.6412	.6155		.6320		.5998
O2	.5346	.2021	.0012	.0015		.1352		.1615
H2	.2871	.0010	.0011	.0011		.0011		.0011
DEF	-.0007		-.0000	.0000	-.0000			
IAT	.0010		.0002	-.0000	-.0108			
UES	-.0003		-.0000	-.0000	.0025			
PCO2	37.9715	37.9715	48.1741	42.7791	48.1742	48.1741	42.7791	
PO2	113.8677	113.8677	36.6162	46.4015	36.6949	36.6162	46.4015	
PH2	61.1608	61.1608	61.1888	61.1888	61.1717	61.1888	61.1888	
(H+)		37.6556	42.3511	39.0214	44.0727	42.9989	40.0941	
PH		7.4242	7.3731	7.4087	7.3558	7.3665	7.3969	
HBO2		.1985				.1340	.1600	
TRANSPORT TIMES								
	AB	UB	UT	AT	AC	UI		
	.1970	.1114	.5912	.3170	.1877	6.0593		
	VE	Q	FB	DERIVATIVES				
	5.9291	6.0000	.7491	.0000	.0003			
PESP FREQ	12.1685	7.5753	1.5811	66.1070	43.9640	DSVOL		
						.1519		

ORIGINAL PAGE IS  
OF POOR QUALITY

EXRT

GRODINS: RESPIRATORY CONTROL MODEL

ADD DATA...  
>@ADD GRODATT

\*RESPIRATORY CHEMOSTAT -- INPUT DATA\*

1	.1783	.5346	.2871	.6412	.0012
6	.0011	.6155	.0015	.0011	6.0000
11	.7491	48.1742	36.6949	61.1717	40.0000
16	.0000	.2000	.1000	.1000	1.1300
21	1.1540	3.0000	1.0000	39.0000	.0500
26	.0500	81.9900	4.3610	2.5240	260.0000
31	.0192	.7000	.2000	.1000	.0000
36	.0078	87.5500	5.3900	.2500	.0000
41	.5470	.5850	.5850	.5850	
45	.1820	.2150	.0000	.0000	

DO YOU WANT GRAPHIC INSTEAD OF TABULAR OUTPUT? (Y/N)  
>Y

INPUT WORK CARDS...  
WORK= WORK LOAD(WATTS)...  
MINS= TIME FOR WORK LOAD...  
PRINT= TIME INCRIMENT(MINS)FOR PRINTOUT...  
EXEC  
MORE= INPUT MORE BEFORE EXEC...  
PUN = EXEC WITH ABOVE THEN CAN INPUT AGAIN...  
STOP= EXEC WITH ABOVE THEN STOP...  
BACK= ERASE PREVIOUS WORK RECORD...  
WORK MINS PRINT EXEC (F6.2,1X,F6.2,1X,F6.2,1X,A4)...  
>@ 2. .5 MORE  
.00 2.00 .50MORE  
WORK MINS PRINT EXEC (F6.2,1X,F6.2,1X,F6.2,1X,A4)...  
>200. 1. .1 MORE  
200.00 1.00 .10MORE  
WORK MINS PRINT EXEC (F6.2,1X,F6.2,1X,F6.2,1X,A4)...  
>200. 7 .5 RUN  
200.00 7.00 .50RUN  
TYPE SHIFT-OUT (SO) AND RETURN-->


ORIGINAL PAGE IS  
OF POOR QUALITY

```

GRAPHIC OUTPUT(Y,N,S),TIME INTERVALS,STARTX,STOPX,(A2,3F5.0)..
>Y 5. 0. 10.
PC02 Y SCALE (A4,8X,F4.0,2F6.0)
PLOT(Y,N,S) LOC HIGH LOW ...
>Y 1. 80. 20.
P02 Y SCALE (A4,8X,F4.0,2F6.0)
PLOT(Y,N,S) LOC HIGH LOW ...
>Y 1. 140. 70.
WORK2 Y SCALE (A4,8X,F4.0,2F6.0)
PLOT(Y,N,S) LOC HIGH LOW ...
>Y 2. 300. 0.
(H+) Y SCALE (A4,8X,F4.0,2F6.0)
PLOT(Y,N,S) LOC HIGH LOW ...
>Y 2. KKK
>Y 3. 54. 42.
UI Y SCALE (A4,8X,F4.0,2F6.0)
PLOT(Y,N,S) LOC HIGH LOW ...
>Y 2. 90. 0.
QB Y SCALE (A4,8X,F4.0,2F6.0)
PLOT(Y,N,S) LOC HIGH LOW ...
>Y 4. 2. 0.
FREQ Y SCALE (A4,8X,F4.0,2F6.0)
PLOT(Y,N,S) LOC HIGH LOW ...
>Y 5. 50. 0.
A0020F Y SCALE (A4,8X,F4.0,2F6.0)
PLOT(Y,N,S) LOC HIGH LOW ...
>Y 6. 200. 0.

```

## GRODINS RESPIRATORY CONTROL MODEL



HARD COPY WANTED(Y/N)...

Y



ORIGINAL PAGE 1\$  
OF POOR QUALITY

>N GRAPHIC OUTPUT(Y,N,S).TIME INTERVALS,STARTX,STOPX,(A2,3F5.0)..

FINAL VALUES FOR FOLLOWING VARIABLES.

1	.20746
2	.52266
3	.26988
4	.66543
5	.00122
6	.00100
7	.77849
8	.00054
9	.00100
10	18.96434
11	.77073
12	48.04239
13	36.28547
14	60.42884

ORIGINAL PAGE IS  
OF POOR QUALITY

APPENDIX C

BATCH EXAMPLE

APPENDIX C  
BATCH EXAMPLE

\*RESPIRATORY CHEMOSTAT -- INPUT DATA\*

1	.1783	.5336	.2881	.6413	.0012
6	.0011	.6153	.0015	.0012	6.0000
11	.7496	48.1202	36.6316	70.6804	40.0000
16	.0000	.2000	.1000	.1000	1.1380
21	1.1540	3.0000	1.0000	39.0000	.0500
26	.0500	81.9900	4.3610	2.5240	260.0000
31	.0192	.7000	.2808	.1000	.0000
36	.0078	87.5500	5.3900	.2500	.0000
41	.5470	.5850	.5850	.5850	
45	.1820	.2150	.0000	.0000	

ORIGINAL PAGE IS  
OF POOR QUALITY

\*\*\*\*\*

WORK LOAD CHG.( .00WATTS FOR .50MINS) AT .0000MINS

TIME	.0000												
	CO2	O2	N2	D E R I V A T I V E S			PCO2	PO2	ALV RQ	.9026	RQ DIFF	-.0271	
									PN2	(H+)	PH	H802	
ALVEOLAR	.1783	.5336	.2881	-.0059	.0102	-.0043	37.9779	113.6568	61.3653				
ARTERIAL	.5653	.2021	.0011				37.9779	113.6568	61.3653	37.6596	7.4241	.1985	
BRAIN	.6413	.0012	.0011	.0006	-.0024	-.0000	48.1939	37.8788	64.1026	42.3632	7.3730		
TISSUE	.6153	.0015	.0012	.0001	-.0003	-.0000	42.7330	47.3485	69.9301	38.9926	7.4090		
CSF				.0090	.1717	-.9675	48.1202	36.6316	70.6804	44.0233	7.3563		
V BRAIN	.6312	.1386	.0011				48.1939	37.8788	64.1026	43.0743	7.3658	.1374	
V TISSUE	.5992	.1632	.0012				42.7330	47.3485	69.9301	40.0956	7.3969	.1617	
TRANSPORT TIMES --		AB	VB	VT	AT	AC **	VI	VE	W	FB	DERIVATIVES		
		.1970	.1114	.5913	.3170	.1877	6.0077	5.9116	6.0000	.7496	.0000	-.0025	
RESP FREQ	10.1710	MINUTE VOLUME		7.5496	DEAD SPACE VENT		1.5900	HEART RATE		66.1070			

TIME	.2500											
	CO2	O2	N2	D E R I V A T I V E S			PCO2	PO2	ALV RQ	.9057	RQ DIFF	-.0302
									PN2	(H+)	PH	H802
ALVEOLAR	.1779	.5349	.2872	.0016	.0009	-.0025	37.8988	113.9349	61.1664			
ARTERIAL	.5649	.2022	.0010				37.8988	113.9349	61.1664	37.6076	7.4247	.1986
BRAIN	.6413	.0012	.0011	-.0002	-.0000	-.0000	48.1942	36.6339	63.6564	42.3634	7.3730	
TISSUE	.6153	.0015	.0012	.0000	-.0000	-.0000	42.7362	46.3936	69.6465	38.9946	7.4090	
CSF				.0087	-.0002	-.9970	48.1224	36.6352	70.4346	44.0254	7.3563	
V BRAIN	.6321	.1352	.0011				48.1942	36.6339	63.6564	43.0120	7.3664	.1341
V TISSUE	.5996	.1615	.0012				42.7362	46.3936	69.6465	40.0670	7.3972	.1600
TRANSPORT TIMES --		AB	VB	VT	AT	AC **	VI	VE	W	FB	DERIVATIVES	
		.1971	.1116	.5912	.3169	.1877	5.9578	5.8639	6.0000	.7471	.0000	-.0010
RESP FREQ	10.1710	MINUTE VOLUME		7.4971	DEAD SPACE VENT		1.5862	HEART RATE		66.1070		

\*\*\*\*\*

WORK LOAD CHG.( 40.00WATTS FOR .25MINS) AT .5000MINS

CHANGE IN METABOLIC RATES MRCO2= .1892 MR02= .2150

TIME	.5000											
	CO2	O2	N2	D E R I V A T I V E S			PCO2	PO2	ALV RQ	.9010	RQ DIFF	-.0255
									PN2	(H+)	PH	H802
ALVEOLAR	.1781	.5352	.2867	-.0008	.0020	-.0013	37.9345	113.9878	61.0777			
ARTERIAL	.5650	.2022	.0010				37.9345	113.9878	61.0777	37.6313	7.4245	.1986
BRAIN	.6412	.0012	.0011	-.0002	.0000	-.0000	48.1800	36.6644	63.2613	42.3547	7.3731	
TISSUE	.6153	.0015	.0012	-.0000	-.0000	-.0000	42.7361	46.2715	69.3700	38.9946	7.4090	
CSF				.0068	.0040	-1.0180	48.1244	36.6357	70.1824	44.0271	7.3563	
V BRAIN	.6320	.1353	.0011				48.1800	36.6644	63.2613	43.0049	7.3665	.1342
V TISSUE	.5997	.1613	.0012				42.7361	46.2715	69.3700	40.0629	7.3973	.1598

APPENDIX C  
BATCH EXAMPLE

TRANSPORT TIMES --	AB	VB	VT	AT	AC **	VI	VE	Q	FB	DERIVATIVES
RESP FREQ 10.1710	.1971	.1116	.5911	.3169	.1877	5.9726	5.8738	6.0000	.7480	.0000 .0005
	MINUTE VOLUME		7.5102	DEAD SPACE	VENT	1.5870	HEART RATE		66.1070	

\*\*\*\*\*  
WORK LOAD CHG.(100.00WATTS FOR .25MINS) AT .7500MINS

CHANGE IN METABOLIC RATES MRCO2= .4367 MR02= .4963

TIME .7500								ALV RQ	1.6623	RQ DIFF	-.7868
CO2	O2	N2	D E R I V A T I V E S			PCO2	P02	PN2	(H+)	PH	HB02
ALVEOLAR .1485	.5788	.2727	.0138	.0227	-.0365	31.6311	123.2812	58.0877			
ARTERIAL .5283	.2033	.0010				31.6311	123.2812	58.0877	33.3900	7.4764	.1994
BRAIN .6413	.0010	.0011	-.0124	-.0003	-.0000	48.1834	32.1767	62.9292	42.3568	7.3731	
TISSUE .6160	.0011	.0012	.0039	-.0010	-.0000	42.8685	35.0427	69.0534	39.0773	7.4081	
CSF			.0068	-.6030	-1.0292	48.1274	36.5570	69.9266	44.0299	7.3563	
V BRAIN .6355	.1217	.0011				48.1834	32.1767	62.9292	42.7587	7.3690	.1207
V TISSUE .6071	.1345	.0012				42.8685	35.0427	69.0534	39.6795	7.4014	.1334
TRANSPORT TIMES --	AB	VB	VT	AT	AC **	VI	VE	Q	FB	DERIVATIVES	
RESP FREQ 12.3689	.1771	.1200	.5517	.2788	.1644	15.8417	16.7335	7.3924	.6117	1.5630	-.1929
	MINUTE VOLUME		18.9704	DEAD SPACE	VENT	2.6828	HEART RATE		78.4258		

TIME 1.0000								ALV RQ	1.3762	RQ DIFF	-.5007
CO2	O2	N2	D E R I V A T I V E S			PCO2	P02	PN2	(H+)	PH	HB02
ALVEOLAR .1536	.5774	.2690	.0108	-.0202	.0094	32.7114	122.9827	57.3059			
ARTERIAL .5350	.2032	.0010				32.7114	122.9827	57.3059	34.1311	7.4668	.1993
BRAIN .6379	.0010	.0011	-.0110	.0000	-.0000	47.4443	31.8706	62.3184	41.9056	7.3777	
TISSUE .6161	.0008	.0012	.0006	-.0005	-.0000	42.8917	26.4801	68.5792	39.0919	7.4079	
CSF			-.0820	-.6229	-1.0804	48.1176	36.3957	69.6638	44.0210	7.3563	
V BRAIN .6323	.1212	.0011				47.4443	31.8706	62.3184	42.2958	7.3737	.1202
V TISSUE .6146	.1052	.0012				42.8917	26.4801	68.5792	39.1903	7.4068	.1044
TRANSPORT TIMES --	AB	VB	VT	AT	AC **	VI	VE	Q	FB	DERIVATIVES	
RESP FREQ 15.2405	.1472	.1195	.4364	.2184	.1345	17.9397	18.6461	9.2192	.6102	2.7336	.0437
	MINUTE VOLUME		21.4438	DEAD SPACE	VENT	3.1508	HEART RATE		94.5195		

ORIGINAL PAGE IS  
OF POOR QUALITY

APPENDIX C  
BATCH EXAMPLE

FINAL VALUES FOR FOLLOWING VARIABLES.

1	.15357
2	.57738
3	.26904
4	.63787
5	.00101
6	.00107
7	.61609
8	.00084
9	.00118
10	9.21921
11	.61019
12	48.11764
13	36.39575
14	69.66384

ORIGINAL PAGE 13  
OF POOR QUALITY

## APPENDIX D

OUTPUT VARIABLES  
NORMAL VALUES

Normal Variables Values are for conditions:

$$B = 260 \text{ mm Hg}$$

$$F_{I(\text{CO}_2)} = .0192$$

$$F_{I(\text{O}_2)} = .7000$$

$$F_{I(\text{N}_2)} = .2808$$

Physiological Variable	Computer Variable	Normal Value	Physiological Variable	Computer Variable	Normal Value
$F_{A(\text{CO}_2)}$	C(1)	.1783	$P_{\text{CSF}(\text{CO}_2)}$	C(12)	48.1743
$F_{A(\text{O}_2)}$	C(2)	.5346	$P_{\text{CSF}(\text{O}_2)}$	C(13)	36.6950
$F_{A(\text{N}_2)}$	C(3)	.2871	$P_{\text{CSF}(\text{N}_2)}$	C(14)	61.1718
$P_a(\text{CO}_2)$	F(7)	37.9705	$C_a(\text{CO}_2)$	CC(1)	.5652
$P_a(\text{O}_2)$	F(1)	113.8681	$C_a(\text{O}_2)$	F(9)	.2021
$P_a(\text{N}_2)$	PAN2	61.1614	$C_a(\text{N}_2)$	F(10)	.0010
$P_B(\text{CO}_2)$	CPB	48.1745	$C_B(\text{CO}_2)$	C(4)	.6412
$P_B(\text{O}_2)$	F(17)	36.6954	$C_B(\text{O}_2)$	C(5)	.0012
$P_B(\text{N}_2)$	F(18)	61.1616	$C_B(\text{N}_2)$	C(6)	.0010
$P_T(\text{CO}_2)$	CPT	42.7786	$C_T(\text{CO}_2)$	C(7)	.6155
$P_T(\text{O}_2)$	PTO2	46.2605	$C_T(\text{O}_2)$	C(8)	.0015
$P_T(\text{N}_2)$	PTN2	61.1673	$C_T(\text{N}_2)$	C(9)	.0010
$C_a(\text{HbO}_2)$	CHB(1)	.1985	$\tau_{aB}$	AB	.1970
$C_{\text{VT}}(\text{HbO}_2)$	CHB(3)	.1342	$\tau_{aT}$	AT	.3170
$C_{\text{VB}}(\text{HbO}_2)$	CHB(2)	.1597	$\tau_{vB}$	VB	.1114

Physiological Variable	Computer Variable	Normal Value	Physiological Variable	Computer Variable	Normal Value
$C_a(H^+)$	CH(1)	37.6549	$\tau_{VT}$	VT	.5912
$C_B(H^+)$	CH(2)	42.3514	$\tau_{ao}$	AC	.1877
$C_T(H^+)$	CH(3)	38.0212	VI	VI	6.0586
$C_{CSF}(H^+)$	CH(4)	44.0728	VE	VE	5.9248
$C_{VB}(H^+)$	HVB	43.0032	Q	C(10)	6.0000
$C_{VT}(H^+)$	HVT	40.0892	FB	C(11)	.7491
$pH_a$	CPH(1)	7.4242	RESP FREQ	FREQ	12.1637
$pH_{CSF}$	PHCSF	7.3558	MINUTE VOL	TVNT	7.5723
$pH_{VB}$	PHVB	7.3665	DEAD SPACE VENT	DEADVT	1.5806
$pH_{VT}$	PHVT	7.3970	HEART RATE	HRATE	66.1070
$pH_{BRAIN}$	CPH(2)	7.3731	ALVEOLAR RQ	RQ	.8754
$pH_{TISSUE}$	CPH(3)	7.4087	RQ DIFF	QF(5)	.0001
AVO2DF	AVO2DF	44.1667	WORK LOAD	WORK2	0.
DSVOL	DSVOL	.1518	MRTCO <sub>2</sub>	RMT(1)	.1820
			MRTO <sub>2</sub>	RMT(2)	.2150



APPENDIX E  
PROGRAM LISTING

## B6-G03432\*TPF5.GRODIN

1		DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
2	1	SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
3	2	RC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
4	3	BQ(4)
5	C	C(40)
6	C	ALVEOLAR VOL GAS FUNCTIONS
7	C	1 FA(CO2)
8	C	2 FA(O2)
9	C	3 FA(N2)
10	C	
11	C	GAS CONCENTRATIONS IN BRAIN.
12	C	4 CB(CO2)
13	C	5 CB(O2)
14	C	6 CB(N2)
15	C	
16	C	GAS CONCENTRATIONS IN TISSUE.
17	C	7 CT(CO2)
18	C	8 CT(O2)
19	C	9 CT(N2)
20	C	CARDIAC OUTPUT.
21	C	10 Q
22	C	CEREBRAL BLOOD FLOW.
23	C	11 QB
24	C	GAS TENSION IN CSF.
25	C	12 PCSF(CO2)
26	C	13 PCSF(O2)
27	C	14 PCSF(N2)
28	C	
29	C	LENGTH OF SIMULATION RUN.
30	C	(THIS IS NOT USED IN TTY MODE. IN BATCH, A WORK CARD WITH 0 TIME WILL
31	C	ALSO STOP RUN).
32	C	15 TMAX
33	C	WEIGHTING OF H+CONC IN CSF VERSUS VENOUS BLOOD OF BRAIN.
34	C	16 CENTRAL SENSITIVITY PARTITION
35	C	BLOOD OXYGEN CAPACITY
36	C	17 (HB)
37	C	TIME CONSTANTS IN CARDIAC OUTPUT AND CEREBRAL BLOOD FLOW RESPONSES.
38	C	18 R1
39	C	19 R2
40	C	
41	C	CONTROLLER EQUATION SENSITIVITY WEIGHTINGS.
42	C	20 CENTRAL SENSITIVITY COEFFICIENT
43	C	21 CAROTID BODY SENSITIVITY COEFFICIENT
44	C	
45	C	VOLUMES OF LUNG, BRAIN, AND TISSUE
46	C	22 KL
47	C	23 KB
48	C	24 KT
49	C	
50	C	BRAIN METABOLIC RATE OF CO2 PRODUCTION.
51	C	25 MRB(CO2)
52	C	BRAIN METABOLIC RATE OF O2 CONSUMPTION.
53	C	26 MRB(O2)
54	C	GAS DIFFUSION COEFF. FOR BLOOD-BRAIN BARRIER.
55	C	27 DCO2
56	C	28 DO2

57	C	29	DN2
58	C		
59	C		BAROMETRIC PRESSURE.
60	C	30	B
61	C		VOL.FRACTION OF INSPIRED GAS.
62	C	31	FI(CO2)
63	C	32	FI(O2)
64	C	33	FI(N2)
65	C		
66	C		VOL.OF CSF.
67	C	34	KCSF
68	C		INITIAL TIME
69	C	35	T
70	C		COMPUTER TIME STEP.
71	C	36	H
72	C		CONTROLLER EQUATION CONSTANT(MAINTAINS RESTING PA(CO2) APPROX.40).
73	C	37	VI(N)
74	C		VALUE FOR RESTING ALVEOLAR VENTILATION.
75	C	38	VI(SS)
76	C		OUTPUT PRINT INCREMENTS (ALSO PRINTS AT .5MIN.INCRIMENTS).
77	C	39	PRINT-ALL TIME
78	C		
79	C		SV(18,50)
80	C		ARTERIAL GAS CONCENTRATIONS AT LUNG EXIT.
81	C	1	CA(CO2)
82	C	2	CA(O2)
83	C	3	CA(N2)
84	C		
85	C		VENOUS GAS CONCENTRATIONS AT BRAIN EXIT.
86	C	4	CVB(CO2)
87	C	5	CVB(O2)
88	C	6	CVB(N2)
89	C		
90	C		VENOUS GAS CONCENTRATIONS AT TISSUE EXIT.
91	C	7	CVT(CO2)
92	C	8	CVT(O2)
93	C	9	CVT(N2)
94	C		
95	C		CARDIAC OUTPUT.
96	C	10	Q
97	C		CEREBRAL BLOOD FLOW.
98	C	11	QB
99	C		TISSUE BLOOD FLOW.
100	C	12	QT
101	C		ARTERIAL H+ CONCENTRATION.
102	C	13	CA(H+)
103	C		ARTERIAL O2 TENSION.
104	C	14	PA(O2)
105	C		
106	C	15	--
107	C		TOTAL GAS CONCENTRATIONS AT BRAIN EXIT.
108	C	16	CVB(CO2) + CVB(O2) + CVB(N2)
109	C		TOTAL GAS CONCENTRATIONS AT TISSUE EXIT.
110	C	17	CVT(CO2) + CVT(O2) + CVT(N2)
111	C		TIME.
112	C	18	T
113	C		

114	C	VTRAN(18)
115	C	ARTERIAL GAS CONCENTRATIONS AT BRAIN ENTRANCE.
116	C	1 CAB(CO2) = CA(CO2)(T = TAB)
117	C	2 CAB(O2) = CA(O2)(T = TAB)
118	C	3 CAB(N2) = CA(N2)(T = TAB)
119	C	
120	C	VENOUS BRAIN GAS CONCENTRATION AT LUNG ENTRANCE.
121	C	4 CVB(CO2)(T = TVB)
122	C	5 CVB(O2)(T = TVB)
123	C	6 CVB(N2)(T = TVB)
124	C	
125	C	VENOUS TISSUE GAS CONCENTRATION AT LUNG ENTRANCE.
126	C	7 CVT(CO2)(T = TVT)
127	C	8 CVT(O2)(T = TVT)
128	C	9 CVT(N2)(T = TVT)
129	C	
130	C	ARTERIAL GAS CONCENTRATIONS AT TISSUE ENTRANCE.
131	C	10 CAT(CO2) = CA(CO2)(T = TAT)
132	C	11 CAT(O2) = CA(O2)(T = TAT)
133	C	12 CAT(N2) = CA(N2)(T = TAT)
134	C	
135	C	ARTERIAL H+ CONCENTRATION AT CAROTID BODIES'SITE.
136	C	13 CAO(H+) = CA(H+)(T = TAO)
137	C	ARTERIAL O2 TENSION AT CAROTID BODIES'SITE.
138	C	14 PAO(O2) = PA(O2)(T = TAO)
139	C	ARTERIAL H+ CONCENTRATION AT BRAIN ENTRANCE.
140	C	15 CAB(H+) = CA(H+)(T = TAB)
141	C	TOTAL GAS CONCENTRATION FROM BRAIN AT LUNG ENTRANCE.
142	C	16 (CVB(CO2) + CVB(O2) + CVB(N2))(T = TVB)
143	C	TOTAL GAS CONCENTRATION FROM TISSUE AT LUNG ENTRANCE.
144	C	17 (CVT(CO2) + CVT(O2) + CVT(N2))(T = TVT)
145	C	
146	C	D(15)
147	C	FOR D(15) THE SYMBOLS B=BAROMETRIC PRESSURE, 47=WATER VAPOR PRESS.,
148	C	K=CONVERSION FACTOR FOR ATM TO MMHG, A=SOLUBILITY COEFF.OF GASES,
149	C	H=COMPUTER TIME STEP, HB=BLOOD OXYGEN CAPACITY
150	C	1 B = 47
151	C	2 K AC02
152	C	3 K AO2
153	C	4 K AN2
154	C	5 K AN2 (B = 47)
155	C	6 K AO2 (B = 47)
156	C	7 K AN2 (B = 47)
157	C	8 0.16 + 2.3(HB)
158	C	9 863/(B = 47)
159	C	10 0.62
160	C	11 K ACSF(CO2)
161	C	12 K ACSF(O2)
162	C	13 K ACSF(N2)
163	C	14 2*H
164	C	15 1.99*H
165	C	F(20)
166	C	COMPARTMENTAL GAS TENSIONS AND CONCENTRATIONS.
167	C	1 PA(O2)
168	C	2 K AC02 PA(CO2)
169	C	3 PB(O2)
170	C	4 K AC02 PB(CO2)

ORIGINAL PAGE IS  
OF POOR QUALITY

```

171      C      5      PT(02)
172      C      6      K AC02 PT(02)
173      C      7      PA(02)
174      C      8      PA(02)
175      C      9      CA(02)
176      C     10      CA(N2)
177      C     11      CA(02) + CA(02) + CA(N2)
178      C     12      CVB(02)
179      C     13      CVT(02)
180      C
181      C  PRODUCT OF DIFFUSION COEFFS. AND GAS DIFFERENTIALS ACROSS BLOOD-BRAIN
182      C  BARRIER.
183      C     14      DC02 (PB(02) - PCSF(02))
184      C     15      DO2 (PB(02) - PCSF(02))
185      C     16      DN2 (PB(N2) - PCSF(N2))
186      C
187      C     17      PB(02)
188      C     18      PB(N2)
189      DIMENSION XNB(4,2), DJ(4), IDJ(2)
190      COMMON/2/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
191      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DT,
192      2      IRK, LOC, ITERX, INDEX, I, J, M, N
193      COMMON/R/ XDS, XMH, CAT, WORK, DUM1, DUM2, DUM3, WORK2, RMTB, RMTB2, TIMEOF
194      1      , RMLIN, ITTY
195      C  ITTY = FLG FOR ITY MODE.
196      C  0 = OUTPUT TO PRINTER (BATCH MODE).
197      C  'TTY' = TTY I/O AND 1ST TIME TO SUBROUTINE RC12.
198      C  1 = TTY I/O AND NOT 1ST TIME TO RC12.
199      DATA ITTY/'TTY'/
200      C  DATA FOR INITIAL CONDITIONS
201      WRITE (6,5)
202      5 FORMAT (' GRODINS: RESPIRATORY CONTROL MODEL'//)
203      300 CONTINUE
204      WRITE(6,483)
205      483 FORMAT('OADD DATA...')
206      C  READ INDICATION OF BATCH OR TTY MODE.
207      READ(5,480) ITTY
208      480 FORMAT(A4)
209      IF(ITTY.NE. ITTTY) ITTY = 0
210      WRITE(6,90)
211      90 FORMAT (1H1,1X,37H*RESPIRATORY CHEMOSTAT *- INPUT DATA*/)
212      C  DATA FOR INITIAL CONDITIONS
213      DO 10 I = 1,40
214      C  1106 HAS PROBLEM WITH END= , SO THIS ISNT USED TO
215      C  DETERMINE END OF RUN(NO CAPABILITY TO START ANOTHER
216      C  MODEL RUN IN SAME COMPUTER RUN).
217      READ(5,190,END=301) C(I), (XN(I,J), J=1,2)
218      10 CONTINUE
219      C  ESTABLISH COMPUTER STEP INDEPENDENT OF INPUT DATA.
220      C(36) = .78125E-2
221      190 FORMAT (5X,F15.0,5X,2A4)
222      DO 20 I = 1,4
223      IP40 = I + 40
224      READ (5,190) BC(I), (XNB(I,J), J = 1,2)
225      20 CONTINUE
226      DO 30 I = 1,2
227      READ (5,190) RMT(I), (XNB(I,J), J = 1,2)

```

```

228      IP40 = I + 44
229      30 CONTINUE
230      DO 40 I = 1,2
231      READ (5,190) DJ(I), (XNR(I,J), J = 1,2)
232      IP40 = I + 46
233      40 CONTINUE
234      C
235      C OUTPUT INPUT DATA.
236      J = 1
237      DO 75 I = 1,8
238      JX = J + 4
239      WRITE(6,92) J, (C(I2), I2=J, JX)
240      92 FORMAT(' ', I2, 2X, 5(F9.4))
241      J = J + 5
242      75 CONTINUE
243      WRITE(6,92) J, (BC(I), I=1,4)
244      J = 45
245      WRITE(6,92) J, RMT(1), RMT(2), DJ(1), DJ(2)
246      C
247      COMMON/PLTRUF/NBUF,XBUF(181),YBUF(181,8),NA(8),KSTOPP,TMAX,KPLT
248      DATA KY/1HY/
249      TMAX=C(15)
250      *WRITE (6,42)
251      42 FORMAT(/'DOO YOU WANT GRAPHIC INSTEAD OF TABULAR OUTPUT? (Y/N)')
252      READ (5,44) KYY
253      44 FORMAT(1A1)
254      IF (KYY.EQ.KY) KPLT=1
255      C IF TTY I/O MAX*TIME WILL COME FROM WORK CARD.
256      IF(ITYY.NE.O) C(15) = 9999999999.
257      C
258      C FI(C02)
259      DUM1=C(31)
260      C FI(O2)
261      DUM2=C(32)
262      C FI(N2)
263      DUM3=C(33)
264      WORK=0.
265      WORK2=0.
266      C METABOLIC RATE OF O2 CONSUMPTION IN TISSUE.
267      RMTB=RMT(2)
268      RMTB2=RMT(2)
269      C
270      TIMEOF=0.
271      XDS=0.
272      XMH=10.*C(36)/0.0078125
273      MMM=0
274      201 CONTINUE
275      XDS=XDS+XMH
276      IF(MMM.EQ.1)XDS=XDS*C(36)
277      MMM=1
278      C(35)=0.
279      C(40)=0.
280      C
281      C INITIAL GUESSES FOR ITERATIVE LOOPS
282      C ARTERIAL CONCENTRATION OF CO2.
283      CC(1) = 0.6
284      C BRAIN CONCENTRATION OF CO2.

```

ORIGINAL PAGE IS  
OF POOR QUALITY

```

285      CC(2) = C(4)
286      C   TISSUE CONCENTRATION OF CO2.
287      CC(3) = C(7)
288      C   BRAIN CO2 TENSION.
289      CPB = 50.0
290      C   TISSUE CO2 TENSION.
291      CPT = 50.0
292      IF(XDS.GT.XMH) GOTO202
293      C   SETS VARIOUS CONSTANTS AND AGGREGATES OF CONSTANTS
294      C   TMAX.
295      C(15) = C(15) + .0001
296      C   PRINT ALL TIME.
297      C   C(39) = C(39) + .0001
298      C   FACTOR OF 1-E-7 MULTIPLYING DIFFUSION COEFFICIENTS.
299      DO 200 I = 27,29
300      C(I) = C(I) * 1.E-7
301      200 CONTINUE
302      202 CONTINUE
303      IRK = 1
304      M = 14
305      N = 5
306      IDJ(1) = 0
307      C   SOLUBILITY COEFFICIENTS.
308      C   A(1) = (ALPHA)CO2, A(2) = (ALPHA)O2, A(3) = (ALPHA)N2,
309      C   A(4) = (ALPHA)CO2, A(5) = (ALPHA)O2, A(6) = (ALPHA)N2
310      A(1) = 0.51
311      A(2) = 0.024
312      A(3) = 0.013
313      A(4) = 0.51
314      A(5) = 0.024
315      A(6) = 0.013
316      C   ATM/MMHG CONVERSION FACTOR.
317      SK = 0.00132
318      C   CARBONIC ACID DISSOCIATION CONSTANT.
319      CADK = 795.0
320      C   VOL(1)-VOL(10) = VOLUMES USED IN CALCULATION OF VARIABLE TIME DELAYS.
321      VOL(1) = 0.015
322      VOL(2) = 1.062
323      VOL(3) = 0.188
324      VOL(4) = 0.06
325      VOL(5) = 0.188
326      VOL(6) = 2.94
327      VOL(7) = 0.735
328      VOL(8) = 1.062
329      VOL(9) = 0.008
330      VOL(10) = 1.062
331      C
332      C (METABOLIC RATE OF CO2 IN BRAIN + TISSUE.) / SAME FOR O2
333      QF(6) = (C(25) + RMT(1))/(C(26) + RMT(2))
334      C   B=47
335      D(1) = C(30) - 47.
336      DO 210 I = 2,4
337      C   PRODUCTS OF CONVERSION FACTORS AND SOLUBILITY COEFFICIENTS.
338      D(I) = SK*A(I-1)
339      D(I+9) = SK*A(I+2)
340      C
341      D(I+3) = D(1)*D(I)

```

ORIGINAL PAGE IS  
OF POOR QUALITY

```

342      210 CONTINUE
343      C FACTOR USED IN ESTABLISHING CAICO2
344          D(8) = 0.16 + 2.3*C(17)
345      C
346          D(9) = 863.0/D(1)
347      C FACTOR USED IN ESTABLISHING CB(CO2).
348          D(10) = 0.62
349      C MANIPULATION OF COMPUTER TIME STEP.
350          D(14) = C(36)*2.0
351          D(15) = D(14) - .01*C(36)
352      C
353          CALL RC3
354          CALL RC4
355          CALL RC5 (CPB, F(4), C(4), BC(2))
356          CALL RC21 (CHB(2), F(3), F(4), C(4), CH(2), CPH(2))
357          CALL RC19 (CPB, CHB(2), CC(2), BC(1), F(4))
358          CALL RC5 (CPT, F(6), C(7), BC(3))
359          CALL RC21 (CHB(3), F(5), F(6), C(7), CH(3), CPH(3))
360          CALL RC19 (CPT, CHB(3), CC(3), BC(1), F(6))
361          CALL RC20
362          CALL RC7
363          CALL RC8
364          CALL RC9
365          CALL RC10
366          CALL RC11
367          CALL RC12
368          GO TO 60
369      50 CALL RC15
370          CALL RC16
371      60 CALL RC13
372          CALL RC12
373      C
374          IF (C(35).GE.XMH) GO TO 201
375      C
376          IF (C(35).GT. C(15)) GOTO 80
377          IF (CX1.GT.C(15)) GOTO 80
378      70 CALL RC14
379          UU = AMOD(C(35), D(14))
380          IF (UU.LT..0001 .OR. UU.GT.D(15)) GOTO 50
381          GOTO 60
382      80 IF (KPLT.LT.1) GO TO 76
383          KSTOPP=1
384          CALL PLOT
385      76 WRITE(6,78)
386      78 FORMAT('1 FINAL VALUES FOR FOLLOWING VARIABLES.')
387          IF (C(37).GT. 1.0E-5) GO TO 250
388      220 CTERM = 0.0
389          IF (VTRAN(14) = 104.0) 230, 240, 240
390      230 CTERM = (23.6E-9)*((104.0 - VTRAN(14))**4.9)
391      240 C(37) = C(20)*(C(16)*VTRAN(15) + (1.0 - C(16))*CH(4))
392          + C(21)*VTRAN(13) + CTERM = V1
393          I = 37
394          WRITE(6,192) I, C(I), (XN(I,J), J = 1,2)
395      250 DO 260 I = 1,14
396          WRITE(6,192) I, C(I), (XN(I,J), J = 1,2)
397      260 CONTINUE
398          WRITE (6,194)

```



```
399      WRITE(6,830)
400      830 FORMAT('ONORMAL TERMINATION')
401      301 CONTINUE
402      STOP
403      C 90 FORMAT (1H148X37H•RESPIRATORY CHEMOSTAT -- INPUT DATA•///)
404      C 92 FORMAT (42X13,10XF10.4,10X2A6)
405      C 190 FORMAT (5XF15.0,5X2A6)
406      192 FORMAT(' ',13,2X,F15.5,2X,2A4)
407      194 FORMAT (1H1)
408      END
```

88-G03432-TPF5-RC3

```

1      SUBROUTINE RC3
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      HC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/27 C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPR, CPT, CADK, X, DT,
8      2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      C6969      FORMAT(1H 7HSUB RC3)
10     C      SETS TIME-DEPENDENT EXPRESSIONS
11     C      TISSUE BLOOD FLOW.
12     QF(1) = C(10) - C(11)
13     C      ARTERIAL O2 TENSION.
14     F(1) = D(1)*C(2)
15     C      ARTERIAL CO2 CONCENTRATION.
16     F(2) = D(5)*C(1)
17     C      BRAIN O2 CONCENTRATION / (CONV*FACTOR*SOLUBILITY COEFF*FOR O2)
18     F(3) = C(5)/D(3)
19     C      (CONV*FACTOR*SOLUBILITY COEFF*FOR CO2) * BRAIN CO2 TENSION.
20     F(4) = D(2)*CPB
21     C      TISSUE O2 CONCENTRATION / (CONV*FACTOR*SOLUBILITY COEFF*FOR O2)
22     F(5) = C(8)/D(3)
23     C      (CONV*FACTOR*SOLUBILITY COEFF*FOR CO2) * TISSUE CO2 TENSION.
24     F(6) = D(2)*CPT
25     C      ARTERIAL CO2 TENSION.
26     F(7) = D(1)*C(1)
27     C      ARTERIAL O2 TENSION.
28     F(8) = D(1)*C(2)
29     RETURN
30     END

```

PRY:5 RC4

ORIGINAL PAGE IS  
OF POOR QUALITY

B6-G03432\*TPFS.RC4

```

1      SUBROUTINE RC4
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/Z/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DT,
8      2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      C      ITERATES FOR CC(1), ARTERIAL CO2 CONCENTRATION
10     C6969      FORMAT(1H 7H SUB RC4)
11     410 CALL RC21 (CHR(1), F(1), F(2), CC(1), CH(1), CPH(1))
12     X = (CC(1) - F(2))/(10.01 * F(1))
13     X = RCF1(X)
14     C SEE EQUATION 3-1, X = CA(CO2)
15     X = BC(1) + 0.375 * (C(17) - CHB(1)) + F(2) - D(8) * (X - 0.14)
16     C CC(1) = CA(CO2)
17     CALL RC6 (CC(1))
18     CC(1) = CC(1) + 2.0 * (X - CC(1)) / 3.0
19     C3000 FORMAT(1H ,5HCC(1),5X,E16.6)
20     IF (ITERX) 420, 410, 420
21     420 RETURN
22     END

```

PRT,S RCS

B6-G03432-TPFS-RC5

```

1      SUBROUTINE RC5 (CP, FB, CCB, BHC)
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/27 C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, GF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DT,
8      2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      C      ITERATES FOR BRAIN AND TISSUE PCO2
10     C6969      FORMAT(1H 7H SUB RC5)
11     510 X = (CCB - FB)/(10.01*CP)
12     X = RCF1(X)
13     C      SEE EQUATION 4.1, X = PB(CO2) .
14     X = (-BHC + CCB + D(10)*(X - 0.14))/D(2)
15     C      CP = PB(CO2) .
16     CALL RC6 (CP)
17     CP = CP + (X - CP)/10.0
18     C      CEREBRAL BLOOD FLOW.
19     FB = D(2)*CP
20     C3000 FORMAT(1H ,4HCP= ,E16.6,4HFB= E16.6,5HCCB= E16.6,5HBHC= E16.6)
21     IF (ITERX) 520, 510, 520
22     520 RETURN
23     END

```

PRT,S RC6

B6-G03432\*TPFS.RC6

```

1      SUBROUTINE RC6 (Y)
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/Z/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DT,
8      2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      C      CHECKS CONVERGENCE OF ITERATIVE PROCEDURES
10     C      RC4 : X=CA(C02), Y=CC(1) .
11     C      RC5 : X=PB(C02), Y=CP .
12     C      RC19 : X=CVB(C02), Y=CVC .
13     C6969      FORMAT(1H 7HSUB RC6)
14     ITERX = 0
15     DIFF = ABS ((X - Y)/Y)
16     IF (DIFF - 1.0E-5) 620, 620, 630
17     620 ITERX = 1
18     630 RETURN
19     END

```

PRT,S RC7

B6-G03432\*YFF5.RC7

```

1      SUBROUTINE RC7
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/Z/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DT,
8      2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      COMMON/R/ XDS,XMH,CAT,WORK,DUM1,DUM2,DUM3,WORK2,RMTB,RMTB2,TIMEOF
10     1      ,RMLIN
11     C6969      FORMAT(1H 7HSUB RC7)
12     C      FILLS SV ARRAY WITH INITIAL CONDITIONS
13     CALL RC16
14     IF(XDS.GT.XMH) GOTO2
15     DO 725 I = 1,17
16     DO 720 J = 2,50
17     SV(I,J) = SV(I,1)
18     720 CONTINUE
19     725 CONTINUE
20     2      CONTINUE
21     DO 730 J = 2,50
22     SV(18,J) = SV(18,J - 1) - D(14)
23     730 CONTINUE
24     C3000 FORMAT(1H ,12H18SV S D(14),6(3X,E16.6)/1H ,6(3X,E16.6)/1H ,7(3X,E1
25     C      C6.6))
26     RETURN
27     END

```

PRT:5 RC8

B6-G03432\*TPFS.RC8

```

1      SUBROUTINE RC8
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/Z/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DI,
8      2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      C      CALCULATES TRANSPORT TIMES
10     C      EQUATIONS 8.10 THRU 8.14
11     C6969      FORMAT(1H 7HSUB RC8)
12      DO 870 I = 1,5
13      DT = C(35) - SV(18,1)
14      ND = 1
15      GO TO (810,812,814,816,818), I
16      810 NC = 11
17      NB = 10
18      GO TO 820
19      812 NC = 10
20      NB = 11
21      GO TO 820
22      814 NC = 10
23      NB = 12
24      GO TO 820
25      816 NC = 12
26      NB = 10
27      QA = QF(1)
28      GO TO 822
29      820 QA = C(NC)
30      822 DO 860 J = 1,2
31      GO TO (834,824), J
32      824 NC = NB
33      ND = K + 1
34      IF (K)      826, 826, 832
35      826 IF (NC - 12)      830, 828, 830
36      828 QA = SV(NC,1) - (SV(NC,1) - QF(1))*DT/(C(35) - SV(18,1))
37      GO TO 834
38      830 QA = SV(NC,1) - (SV(NC,1) - C(NC))*DT/(C(35) - SV(18,1))
39      GO TO 834
40      832 QA = SV(NC,ND) - (SV(NC,K) - SV(NC,ND))*DT/D(14)
41      834 IJ = 2*I + J - 2
42      AB = VOL(IJ)
43      AA = DT*(QA + SV(NC,ND))/2.0
44      DO 838 K = ND,49
45      IF (AA - AB)      836, 836, 840
46      836 AA = AA + C(36)*(SV(NC,K) + SV(NC,K+1))
47      838 CONTINUE
48      WRITE (6,890) I
49      840 DA = AA - AB
50      K = K - 1
51      IF (K)      842, 842, 846
52      842 DV = SV(NC,1) - QA
53      IF (DV)      850, 844, 850
54      844 DT = DA/QA
55      GO TO 860
56      846 DV = SV(NC,K+1) - SV(NC,K)

```

ORIGINAL PAGE IS  
OF POOR QUALITY

```
57      IF (DV)      850, 848, 850
58      848 DT = DA/SV(NC,K)
59      GO TO 860
60      850 DT = (SV(NC,K+1) - SQRT (SV(NC,K+1)**2 - DV*DA/C(36)))/(DV/D(14))
61      860 CONTINUE
62      TAU(1) = C(35) - SV(18,K + 1) - DT
63      870 CONTINUE
64      RETURN
65      890 FORMAT (5X27HSV ARRAY EXCEEDED ON CYCLE 12)
66      END
```

PRT,S RC9



ORIGINAL PAGE IS  
OF POOR QUALITY

B6-G03432\*TPFS.RC9

```

1      SUBROUTINE RC9
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DW(4)
6      COMMON/2/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DW, VE, VI, CPB, CPT, CADK, X, DT,
8      2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      C      SETS VALUES IN VTRAN ARRAY
10     C6969      FORMAT(1H 7HSUB RC9)
11     DO 960 I = 1,5
12     TA = TAU(I) = (C(35) - SV(18,1))
13     LOC = TA/D(14)
14     IF (LOC = 49)      904, 904, 902
15     902 WRITE (6,990) I,LOC
16     LOC = 49
17     904 XLLOC = LOC
18     TB = XLLOC*D(14)
19     DT = TA = TB
20     GO TO (910,920,930,940,950), I
21     910 DO 914 J = 1,3
22     C LUNG TO BRAIN CO2,O2,N2 TIME DELAYED ARTERIAL CONCENTRATIONS.
23     VTRAN(J) = RCF3(J)
24     914 CONTINUE
25     C LUNG TO BRAIN H+ TIME DELAYED ARTERIAL CONCENTRATION.
26     VTRAN(15) = RCF3(13)
27     GO TO 960
28     920 DO 924 J = 4,6
29     C BRAIN TO LUNG CO2,O2,N2 TIME DELAYED VENOUS CONCENTRATIONS.
30     VTRAN(J) = RCF3(J)
31     924 CONTINUE
32     C BRAIN TO LUNG COMBINED CO2,O2,N2 TIME DELAYED VENOUS CONCENTRATIONS.
33     VTRAN(16) = RCF3(16)
34     GO TO 960
35     930 DO 934 J = 7,9
36     C TISSUE TO LUNG CO2,O2,N2 TIME DELAYED VENOUS CONCENTRATIONS.
37     VTRAN(J) = RCF3(J)
38     934 CONTINUE
39     C TISSUE TO LUNG COMBINED CO2,O2,N2 TIME DELAYED VENOUS CONCENTRATIONS.
40     VTRAN(17) = RCF3(17)
41     GO TO 960
42     940 DO 944 J = 1,3
43     C LUNG TO TISSUE CO2,O2,N2 TIME DELAYED ARTERIAL CONCENTRATIONS.
44     VTRAN(J+9) = RCF3(J)
45     944 CONTINUE
46     GO TO 960
47     C LUNG TO CAROTID SITE H+ TIME DELAYED ARTERIAL CONCENTRATION.
48     950 VTRAN(13) = RCF3(13)
49     C LUNG TO CAROTID SITE O2 TIME DELAYED ARTERIAL TENSION.
50     VTRAN(14) = RCF3(14)
51     960 CONTINUE
52     C      NAMELIST/DONM/VTRAN
53     RETURN
54     990 FORMAT (5X27HSV ARRAY EXCEEDED ON CYCLE 12,12H WITH LOC = 14)
55     END

```

R6-G03432\*TPFS\*RC10

ORIGINAL PAGE IS  
OF POOR QUALITY

```

1      SUBROUTINE RC10
2      DIMENSION C(40), XN(40,2), SV(18), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/Z/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DT,
8      2      IRK, LOC, ITERX, INDEX, I, J, K, N
9      C6969      FORMAT(1H 8HSUB RC10)
10     C      COMPUTES EMPIRICAL FUNCTIONS FOR ACRDIAC OUTPUT AND BRAIN BLOOD
11     C      FLOW DIFFERENTIAL EQUATIONS
12     C      F(8) > PA(02) .
13     IF (F(8) - 104.0)      1008, 1020, 1020
14     C      (DELTA)Q(02) , EQUATION 7.3 .
15     1008 DQ(1) = ((-1.0033E-5*F(8) + 2.9241E-3)*F(8) - 0.2885)*F(8) + 9.6651
16     C      (DELTA)QB(02) , EQUATION 7.9 .
17     DQ(2) = (((7.6559E-8*F(8) - 2.324E-5)*F(8) + 2.6032E-3)*F(8)
18     1      - 0.1323)*F(8) + 2.785
19     IF (DQ(1))      1012, 1016, 1016
20     1012 DQ(1) = 0.0
21     1016 IF (DQ(2))      1024, 1028, 1028
22     1020 DQ(1) = 0.0
23     1024 DQ(2) = 0.0
24     C      F(7) = PA(02) .
25     1028 IF (F(7) - 60.0)      1032, 1032, 1036
26     C
27     C      IF PCO2 GT 60 DQ(3) STAYS AT ITS VALUE AT 60 - - OLD ROUTINE SETS
28     C      THE VALUE OF DQ(3) EQUAL TO 0
29     1032 IF (F(7) - 40.0)      2036, 1040, 1040
30     2036 DQ(3)=0.
31     GOTO 1044
32     C      (DELTA)Q(CO2) , REPLACES EQUATION 7.6 .
33     1036 DQ(3)=6.0
34     C
35     GO TO 1044
36     C      (DELTA)Q(CO2) , EQUATION 7.5 .
37     1040 DQ(3) = 0.3*(F(7) - 40.0)
38     1044 IF (F(7) - 38.0)      1048, 1052, 1052
39     C      (DELTA)QB(CO2) , EQUATION 7.11 .
40     1048 DQ(4) = (8.0163E-4*F(7) - 3.1073E-2)*F(7) + 2.3232E-2
41     RETURN
42     1052 IF (F(7) - 44.0)      1056, 1056, 1060
43     1056 DQ(4) = 0.0
44     RETURN
45     C      (DELTA)QB(CO2) , EQUATION 7.13 .
46     1060 DQ(4) = (((-2.1748E-7*F(7) + 9.3918E-5)*F(7) - 1.2947E-2)*F(7)
47     1      + 0.7607)*F(7) - 15.58
48     C      NAMELIST/DG/DQ,F
49     RETURN
50     END

```

PRT:5 RC11

ORIGINAL PAGE 18  
OF POOR QUALITY

186-G03432\*TPF\*.RC11

```

1      SUBROUTINE RC11
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/27 C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DT,
8      2      IRK, LOC, ITERX, INDEX, I, J, N, N
9      C      CALCULATES DIFFERENTIAL EQUATIONS
10     C6969      FORMAT(1H 8H SUB RC11)
11     CALL RC17
12     C EQUATION 10.1 .
13     DC(1) = (VI*C(31) - VE*C(1) + D(9)*(C(11)*VTRAN(4) + QF(1)
14     1      *VTRAN(7) - C(10)*CC(1)))/C(22)
15     C EQUATION 10.2 .
16     DC(2) = (VI*C(32) - VE*C(2) + D(9)*(C(11)*VTRAN(5) + QF(1)
17     1      *VTRAN(8) - C(10)*F(9)))/C(22)
18     C EQUATION 10.3 .
19     DC(3) = (VI*C(33) - VE*C(3) + D(9)*(C(11)*VTRAN(6) + QF(1)
20     1      *VTRAN(9) - C(10)*F(10)))/C(22)
21     C EQUATION 10.4 .
22     DC(4) = (C(25) + C(11)*(VTRAN(1) - CC(2)) - F(14))/C(23)
23     C EQUATION 10.5 .
24     DC(5) = (-C(26) + C(11)*(VTRAN(2) - F(12)) - F(15))/C(23)
25     C EQUATION 10.6 .
26     DC(6) = (C(11)*(VTRAN(3) - C(6)) - F(16))/C(23)
27     C EQUATION 10.7 .
28     DC(7) = (RMT(1) + QF(1)*(VTRAN(10) - CC(3)))/C(24)
29     C EQUATION 10.8 .
30     DC(8) = (-RMT(2) + QF(1)*(VTRAN(11) - F(13)))/C(24)
31     C EQUATION 10.9 .
32     DC(9) = QF(1)*(VTRAN(12) - C(9))/C(24)
33     C EQUATION 7.1 .
34     DC(10) = (-C(10) + 6.0 + DQ(1) + DQ(3))/C(18)
35     C DEPENDANCE OF CARDIAC OUTPUT ON TISSUE
36     C UTILIZATION OF OXYGEN.
37     C
38     XAB=5.5 *(RMT(2)-.215)+6.-C(10)
39     IF((RMT(2).GT..215).AND.(XAB.GT.0.))DC(10)=DC(10)+XAB/.010
40     C
41     C
42     C EQUATION 7.7 .
43     DC(11) = (-C(11) + 0.75 + DQ(2) + DQ(4))/C(19)
44     C EQUATION 1.10 .
45     DC(12) = F(14)/(C(34)*D(11))
46     C EQUATION 1.11 .
47     DC(13) = F(15)/(C(34)*D(12))
48     C EQUATION 1.12 .
49     DC(14) = F(16)/(C(34)*D(13))
50     C NAMELIST/AB/DC
51     RETURN
52     END

```

PRT:5 RC12

ORIGINAL PAGE IS  
OF POOR QUALITY

B6-G03432\*TPF5.RC12

```

1      SUBROUTINE RC12
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/2/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DT,
8      2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      COMMON/R/ XDS,XMH,CXT,WORK,DUM1,DUM2,DUM3,WORK2,RMTB,RMTB2,TIMEOF
10     1      ,RMLIN,ITTY
11     COMMON/PLTBUF/NBUF,XBUF(181),YBUF(181,8),NA(8),KSTOPP,TMAX,KPLT
12     DATA NA/' PC02  P02 WORK2 (H+)  VI  QB  FREQAV02DF'/
13     DATA IRUN/'RUN '/,ISTOP/'STOP'/,MORE/'MORE'/
14     DATA IBACK/'BACK'/
15     DIMENSION WRKTTY(50,3)
16     C6969      FORMAT(1H 8H SUB RC12)
17     C      OUTPUT -- PUNCHED CARDS AND PRINTED
18     CXT=C(35)+XDS-10.
19     IF(CXT.LE.0.)CXT=+0.
20     C      DEAD SPACE VOLUME
21     DSVOL=0.140+0.002*VE
22     C      RESPIRATORY FREQUENCY.
23     FREQ=((1.+(.726*VE)/DSVOL)**.5-1.)/.363
24     C      DEAD SPACE VENTILATION
25     DEADVT=1.+.098*VE
26     C      C(31)=(DEADVT*C(1)+VE*DUM1)/(DEADVT+VE)
27     C      C(32)=(DEADVT*C(2)+VE*DUM2)/(DEADVT+VE)
28     C      C(33)=(DEADVT*C(3)+VE*DUM3)/(DEADVT+VE)
29     C      MINUTE VOLUME.
30     TVNT=DEADVT*(VE+VI)/2.
31     C      HEART RATE.
32     HRATE=43.8*(RMT(2)+C(26))+54.5
33     C
34     C
35     IF(CXT .LT. TIMEOF) GO TO 203
36     C
37     C      HERE IF NEED TO READ A NEW WORK LOAD CARD.
38     C      BRANCH IF IN BATCH MODE.
39     IF(ITTY .EQ. 0) GO TO 500
40     C
41     C
42     C      HERE IF TTY MODE.....
43     IF(ITTY .EQ. 1) GO TO 550
44     C      HERE IF TTY MODE, AND 1ST TIME THIS ROUTINE CALLED.
45     ITTY = 1
46     WRITE(6,505)
47     505 FORMAT('OINPUT WORK CARDS...')
48     1 ' WORK= WORK LOAD(WATTS)...'
49     2 ' MINS= TIME FOR WORK LOAD...'
50     3 ' PRINT= TIME INCRIMENT(MINS)FOR PRINTOUT...'/
51     4 ' EXEC...'/
52     5 ' MORE= INPUT MORE BEFORE EXEC...'/
53     6 ' RUN = EXEC.WITH ABOVE,THEN CAN INPUT AGAIN...'/
54     7 ' STOP= EXEC.WITH ABOVE THEN STOP...'/
55     8 ' BACK= ERASE PREVIOUS WORK RECORD...')
56     504 ITTYIN = 0

```

```

57      ITTYOT = J
58      C
59      501 IF(ITTYIN .LT. 50) GO TO 506
60      C HERE IF BUFFER FOR WORK LOAD CARDS IS FULL.
61      WRITE(6,511)
62      511 FORMAT('0BUFFER FOR WORK LOAD RECORDS FULL.')
63      I ' WILL USE EXEC= RUN.')
64      LEXEC = IRUN
65      GO TO 551
66      C
67      506 ITTYIN = ITTYIN + 1
68      IF (KSKP.GT.1 .AND. KPLT.GT.0) CALL PAGE3
69      509 WRITE(6,507)
70      507 FORMAT(' WORK MINS PRINT EXEC ',
71      I '(F6.2,1X,F6.2,1X,F6.2,1X,A4),.,.')
72      READ(5,502,ERR=509) (WRKTTY(ITTYIN,J),J=1,3),LEXEC
73      502 FORMAT(F6.2,1X,F6.2,1X,F6.2,1X,A4)
74      WRITE(6,503) (WRKTTY(ITTYIN,J),J=1,3),LEXEC
75      503 FORMAT(3(' ',F6.2),A4)
76      IF(LEXEC .NE. 1BACK) GO TO 518
77      ITTYIN = ITTYIN - 1
78      IF(ITTYIN .LT. 1) ITTYIN = 1
79      GO TO 509
80      C
81      518 IF(LEXEC .EQ. IRUN .OR. LEXEC .EQ. 1STOP) GO TO 551
82      IF(LEXEC .EQ. MORE) GO TO 501
83      WRITE(6,510)
84      510 FORMAT(' EXEC.PARAMETER WRONG. TRY AGAIN.')
85      GO TO 509
86      C
87      C HERE IF 1ST TIME THIS ROUTINE CALLED.
88      C SEE IF MORE WORK CARDS IN BUFFER(WRKTTY(50,3))
89      550 IF(ITTYOT .LE. ITTYIN) GO TO 551
90      C HERE IF EXHAUSTED WORK CARD BUFFER (WRKTTY(50,3)).
91      IF(LEXEC .EQ. IRUN) GO TO 504
92      C FORCE END OF COMPUTER RUN WHEN LEXEC= 'STOP'.
93      C(15) = 0.
94      GO TO 1210
95      C
96      551 WORK2 = WRKTTY(ITTYOT,1)
97      DURAT = WRKTTY(ITTYOT,2)
98      C(39) = WRKTTY(ITTYOT,3)
99      ITTYOT = ITTYOT + 1
100     KSTOPP=2
101     IF (KSKP.LT.1 .OR. KPLT.LT.1) GO TO 557
102     CALL PLOT
103     557 KSKP=2
104     KSTOPP=0
105     GO TO 606
106     C
107     C
108     C
109     203 IF(MARKER.EQ.0) GOTO101
110     I WORK=WORK2
111     MARKER=1
112     C SYSTEM RESPONSES: TIME CONSTANTS FOR WORK LOAD LEVELS(INCREASING).
113     IF(WORK.LE.0.)GOTO2

```

```
114      IF(WORK.GE.50.) TCT=2.3/(2.*WORK/200.)
115      IF(WORK.LT.50) TCT=4.6
116      C  TISSUE O2 METABOLIC RATE.
117      RMT(2)=SSO2W(WORK)-(SSO2W(WORK)-RMTB2)*EXP(-TCT*(CXT -TIMEON))
118      VTIME=1.1-1.1*EXP(-TCT*(CXT-TIMEON)/1.92)
119      C  TERM USED IN VI THAT IS A COMPONENT OF TRANSIENT RESPONSE RELATED
120      C  TO WORK LOAD.
121      RMLIN =SSO2W(WORK)-(SSO2W(WORK)-RMTB2)*(1.-VTIME)
122      IF(VTIME.GE.1.) RMLIN=SSO2W(WORK)
123      C  TISSUE CO2 METABOLIC RATE.
124      RMT(1)=.88*RMT(2)
125      IF(TVNT.GT.37.) RMT(1)=(TVNT+40.77)*RMT(2)/88.5
126      IF(C(35).LT.C(40)) GOTO2
127      IF (KPLT.GT.0) GO TO 2
128      WRITE (6,333) RMT(1),RMT(2)
129      333  FORMAT( '0',1X,25HCHANGE IN METABOLIC RATES,5X,7HMRCO2= ,F10.4,
130      1  5X,6HMR02= ,F10.4,/)
131      C
132      C
133      2  CONTINUE
134      C  U = AMOD(C(35), 0.5)
135      C  IF (U .LT. 1.0E-5 .OR. U .GT. .4999)      GO TO 1210
136      IF(C(35).LT.C(40))GOTO1230
137      C(40)=C(40)+C(39)
138      C  ARTERIAL N2 TENSION.
139      1210 PAN2 = D(1)*C(3)
140      C  TISSUE O2 TENSION.
141      PTO2 = C(8)/D(3)
142      C  TISSUE N2 TENSION.
143      PTN2 = C(9)/D(4)
144      C  CEREBROSPINAL FLUID PH , EQUATION 6.2 .
145      PHCSF = 9. - RCF1(CH(4))
146      C  VENOUS BRAIN H+ CONCENTRATION , EQUATION 4.7 .
147      HVB = CADK*F(4)/(CC(2) - F(4))
148      C  VENOUS BRAIN PH , EQUATION 4.6 .
149      PHVB = 9. - RCF1(HVB)
150      C  VENOUS TISSUE H+ CONCENTRATION , EQUATION 5.7 .
151      HVT = CADK*F(6)/(CC(3) - F(6))
152      C  VENOUS TISSUE PH , EQUATION 5.6 .
153      PHVT = 9. - RCF1(HVT)
154      C  RESPIRATORY QUOTIENT (ALVEOLAR).
155      RQ = ((C(11)*VTRAN(4) + QF(1)*VTRAN(7))/C(10) - CC(1))/
156      1  (F(9) - (C(11)*VTRAN(5) + QF(1)*VTRAN(8))/C(10))
157      QF(5) = QF(4) - RQ
158      AVO2DM=(F(9)*C(10)-F(13)*(C(10)-C(11))-F(12)*C(11))*1000.
159      AVO2DF=AVO2DM/C(10)
160      C
161      C
162      C  HERE WHEN READY TO PRINT.
163      C  SEE IF TTY MODE.
164      IF(ITTY .EQ. 0) GO TO 610
165      C
166      NBUF=NBUF+1
167      IF (NBUF.GT.181) NBUF=181
168      XBUF(NBUF)=CXT
169      YBUF(NBUF,1)=F(7)
170      YBUF(NBUF,2)=F(1)
```

ORIGINAL PAGE IS  
OF POOR QUALITY

```

171      YBUF(NBUF,3)=WORK2
172      YRUF(NBUF,4)=CH(4)
173      YBUF(NBUF,5)=VI
174      YRUF(NBUF,6)=C(11)
175      YRUF(NBUF,7)=PREG
176      YRUF(NBUF,8)=AVO2DF
177      IF (KPLT.GT.0) CALL PLOT
178      IF (KPLT.GT.0) GO TO 1230
179      C
180      C  HERE IF TTY OUTPUT.
181      WRITE(6,700) CXT,RQ,QF(5)
182      700 FORMAT(' QTIME',F10.4,' MINS',3X,'ALV RQ',F10.4,
183      1 3X,'RQ DIFF',F8.4/
184      2 7X,'ALVEOLAR ARTERIAL    BRAIN    TISSUE',6X,
185      3 'CSF    V BRAIN V TISSUE')
186      WRITE(6,701) C(1),CC(1),C(4),C(7),CC(2),CC(3)
187      701 FORMAT('  CO2 ',4(F9.4),9X,2(F9.4))
188      WRITE(6,702) C(2),F(9),C(5),C(8),F(12),F(13)
189      702 FORMAT('    O2 ',4(F9.4),9X,2(F9.4))
190      WRITE(6,703) C(3),F(10),C(6),C(9),C(6),C(9)
191      703 FORMAT('    N2 ',4(F9.4),9X,2(F9.4))
192      WRITE(6,704) DC(1),DC(4),DC(7),DC(12)
193      704 FORMAT('  DER  ',F9.4,9X,3(F9.4))
194      WRITE(6,705) DC(2),DC(5),DC(8),DC(13)
195      705 FORMAT('  IVAT ',F9.4,9X,3(F9.4))
196      WRITE(6,706) DC(3),DC(6),DC(9),DC(14)
197      706 FORMAT('  IVES ',F9.4,9X,3(F9.4))
198      WRITE(6,707) F(7),F(7),CPB,CPT,C(12),CPB,CPT
199      707 FORMAT(' PCO2 ',7(F9.4))
200      WRITE(6,708) F(1),F(1),F(17),PTO2,C(13),F(17),PTO2
201      708 FORMAT('  PO2 ',7(F9.4))
202      WRITE(6,709) PAN2,PAN2,F(18),PTN2,C(14),F(18),PTN2
203      709 FORMAT('  PN2 ',7(F9.4))
204      WRITE(6,710) (CH(I),I=1,4),HVB,HVT
205      710 FORMAT(' (H*) ',9X,6(F9.4))
206      WRITE(6,711) (CPH(I),I=1,3),PHCSF,PHVB,PHVT
207      711 FORMAT('  PH ',9X,6(F9.4))
208      WRITE(6,712) (CHB(I),I=1,3)
209      712 FORMAT(' HB02 ',9X,F9.4,27X,2(F9.4))
210      C  PRINTOUT TRANSPORT TIMES.
211      WRITE(6,713) (TAU(I),I=1,5),VI
212      713 FORMAT(' TRANSPORT TIMES'/
213      1 ' ',12X,'AB',7X,'VB',7X,'VT',7X,'AT',7X,'AC',7X,'VI'/
214      2 ' ',5X,6(F9.4))
215      WRITE(6,714) VE,C(10),C(11),DC(10),DC(11)
216      714 FORMAT(' ',12X,'VE',8X,'Q',7X,'FB',4X,'DERIVATIVES'/
217      1 ' ',5X,5(F9.4))
218      WRITE(6,715) FREQ,TVNT,DEADVT,HRATE,AVO2DF,DSVOL
219      715 FORMAT(' ',5X,'RESP FREQ  MIN VOL D.S. VENT  HEART R  AVO2DF  '
220      1 ',  DSVOL'/ 6X,6(F9.4))
221      C
222      RETURN
223      610 IF (N.NE.4) GO TO 1220
224      N = 0
225      WRITE (6,1805)
226      1220 N = N + 1
227      C

```

```

228      WRITE (6,1810) CXT, RQ, QF(5)
229      C
230      WRITE (6,1815) (C(I), I = 1,3), (DC(I), I = 1,3), F(7), F(1),
231      1      PAN2
232      WRITE (6,1820) CC(1), F(9), F(10), F(7), F(1), PAN2, CH(1),
233      1      CPH(1), CHB(1)
234      WRITE (6,1825) (C(I), I = 4,6), (DC(I), I = 4,6), CPB, F(17),
235      1      F(18), CH(2), CPH(2)
236      WRITE (6,1830) (C(I), I = 7,9), (DC(I), I = 7,9), CPT, PT02,
237      1      PTN2, CH(3), CPH(3)
238      WRITE (6,1835) (DC(I), I = 12,14), (C(I), I = 12,14), CH(4),
239      1      PHCSF
240      WRITE (6,1840) CC(2), F(12), C(6), CPB, F(17), F(18), HVB,
241      1      PHVB, CHB(2)
242      WRITE (6,1845) CC(3), F(13), C(9), CPT, PT02, PTN2, HVT,
243      1      PHVT, CHB(3)
244      WRITE (6,1850) (TAU(I), I = 1,5), VI, VE, C(10), C(11), DC(10),
245      1      DC(11)
246      WRITE (6,1855) FREQ,TVNT,DEADVT,HRATE,AV02DF,DSVOL
247      1230 RETURN
248      1290 FORMAT (5H XXXX5X7F10.4)
249      1292 FORMAT (8F10.4)
250      1805 FORMAT (1H1)
251      1810 FORMAT (1H06X4HTIMEF10.4,74X6HALV RQF10.4,3X7HRQ DIFF,F8.4/
252      1      16X3HCO28X2H028X2HN27X21HD E R I V A T I V E S9X4HPC026X
253      2      3HP027X3HPN27X4H(H+)7X2HPM5X4HHB02)
254      1815 FORMAT (3X8HALVEOLAR9F10.4)
255      1820 FORMAT (3X8HARTERIAL3F10.4,30X,5F10.4,F8.4)
256      1825 FORMAT (6X5HBRAIN11F10.4)
257      1830 FORMAT (5X6HTISSUE11F10.4)
258      1835 FORMAT (8X3HCSF30X8F10.4)
259      1840 FORMAT (4X7HV BRAIN3F10.4,30X,5F10.4,F8.4)
260      1845 FORMAT (3X8HV TISSUE3F10.4,30X,5F10.4,F8.4)
261      1850 FORMAT (5X18HTRANSPORT TIMES --4X2HAB8X2HVB8X2HVT8X2HAI8X2HAC2X
262      1      2H**4X2HV18X2HVE8X1HQ9X2HFB7X11HDERIVATIVES/21X,10F10.4,F8.4)
263      1855 FORMAT(3X,9HRESP FREQ,F8.4,2X,13HMINUTE VOLUME,F8.4,
264      1      2X,8HD S VENT,F8.4,2X,10HHEART RATE,F8.4,
265      2      2X,7HAV02DF,F8.4,2X,5HDSVOL,F8.4)
266      C BATCH MODE WORK CARD READ...
267      C
268      C WILL USE WORK CARD WITH TIME=0 AS INDICATION
269      C OF END OF RUN BECAUSE 1106 HAS PROBLEM
270      C WITH END= ON READ.
271      500 READ(5,300,END=2) WORK2,DURAT
272      300 FORMAT(F6.2,3X,F6.2)
273      C
274      IF(DURAT .GT. 0.) GO TO 606
275      C HERE IF READ INDICATION OF END OF RUN IN BATCH MODE.
276      C(15) = 0.
277      GO TO 1210
278      C
279      606 IF (KPLT.GT.0) GO TO 607
280      WRITE (6,305) WORK2,DURAT,CXT
281      305 FORMAT('D',43(' '),/
282      1      ' WORK LOAD CHG.(' ,F6.2,' WATTS FOR',
283      2      F6.2,' MINS) AT',F9.4,' MINS')
284      607 TIMEOF=DURAT+CXT

```



```
285      TIMEON=CXT
286      C  SYSTEM RESPONSES: TIME CONSTANTS FOR WORK LOADS AND TISSUE O2
287      C  METABOLIC RATE.
288          IF(WORK2.GE.WORK)RMTB2=RMT(2)
289      C  DECREASING WORK LOADS.
290          IF(WORK2.LT.WORK) RMTM=RMT(2)
291          IF(WORK2.LT.WORK)RMTB=5502W(WORK2)
292          IF((WORK2.LT.WORK).AND.(WORK.GE.50.)) TCT=2.3/(2.*WORK/200.)
293          IF((WORK2.LT.WORK).AND.(WORK.LT.50.))TCT=4.6
294          IF(WORK2.GE.WORK) GOTO1
295      101 WORK=WORK2
296          MARKER=0
297      C  TISSUE O2 METABOLIC RATE.
298          RMT(2)=RMTB-(RMTB-RMTM)*EXP(-TCT*(CXT-TIMEON)*.50)
299          VTIME=1.-1.-1*EXP(-TCT*(CXT-TIMEON)/3.84)
300      C  TERM USED IN V1 THAT IS A COMPONENT OF TRANSIENT RESPONSE RELATED
301      C  TO WORK LOAD.
302          RMLIN =RMTB-(RMTB-RMTM)*((1.-VTIME)
303          IF(VTIME.GE.1.) RMLIN=RMTB
304      C  TISSUE CO2 METABOLIC RATE.
305          RMT(1)=.88*RMT(2)
306          IF(TVNT.GT.37.) RMT(1)=(TVNT+40.77)*RMT(2)/88.5
307          IF(C(35).LT.C(40)) GOTO2
308          IF (KPLT.GT.0) GO TO 2
309          WRITE (6,333) RMT(1),RMT(2)
310          GOTO2
311      END
```

PRINT,5 RC13

IB6-G03432\*TPFS.RC13

```

1      SUBROUTINE RC13
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3          1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4          2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5          3      DQ(4)
6      COMMON/Z/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7          1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DT,
8          2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      C6969      FORMAT(1H 8HSUB RC13)
10     C      SOLVES M DIFFERENTIAL EQUATIONS BY FOURTH-ORDER RUNGE-KUTTA AND
11     C      ADAMS-MOULTON PREDICTOR-CORRECTOR METHODS
12     C      NAMELIST/DBG/C,DC,SC
13     IF (IRK = 4)      1304, 1356, 1356
14     1304 DO 1352 INDEX = 1,4
15     DO 1308 I = 1,M
16     RK(I,INDEX) = DC(I)
17     1308 CONTINUE
18     GO TO (1312, 1320, 1328, 1340), INDEX
19     1312 DO 1316 I = 1,M
20     SC(I,IRK+1) = C(I)
21     SC(I,IRK) = DC(I)
22     1316 CONTINUE
23     TI = C(35)
24     1320 C(35) = TI + C(36)/2.0
25     DO 1324 I = 1,M
26     C(I) = SC(I,IRK+1) + C(36)*RK(I,INDEX)/2.0
27     1324 CONTINUE
28     GO TO 1336
29     1328 C(35) = TI + C(36)
30     DO 1332 I = 1,M
31     C(I) = SC(I,IRK+1) + C(36)*RK(I,INDEX)
32     1332 CONTINUE
33     1336 CALL RC14
34     GO TO 1352
35     1340 DO 1344 I = 1,M
36     C(I) = SC(I,IRK+1) + C(36)*(RK(I,1) + 2.0*RK(I,2) + 2.0*RK(I,3)
37     I      + RK(I,4))/6.0
38     1344 CONTINUE
39     IRK = IRK + 1
40     1352 CONTINUE
41     RETURN
42     1356 DO 1360 I = 1,M
43     SC(I,5) = C(I)
44     SC(I,4) = DC(I)
45     C(I) = SC(I,5) + C(36)*(55.0*SC(I,4) - 59.0*SC(I,3) + 37.0*SC(I,2)
46     I      - 9.0*SC(I,1))/24.0
47     1360 CONTINUE
48     C(35) = C(35) + C(36)
49     NC35=C(35)/C(36) + .1
50     C(35)=C(36)*NC35
51     1364 CALL RC14
52     DO 1368 I = 1,M
53     SC(I,1) = C(I)
54     C(I) = SC(I,5) + C(36)*(9.0*DC(I) + 19.0*SC(I,4) - 5.0*SC(I,3)
55     I      + SC(I,2))/24.0
56     1368 CONTINUE

```

```
57      DO 1372 I = 1,M
58      IF (ABS (C(I) - SC(I,1)) - 1.0E-3) 1372, 1372, 1364
59      1372 CONTINUE
60      DO 1376 I = 1,M
61      DO 1376 J = 1,3
62      SC(I,J) = SC(I,J+1)
63      1376 CONTINUE
64      RETURN
65      END
```

PRINT, S RC14

DB6-G03432\*TPF5.RC14

```

1      SUBROUTINE RC14
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/2/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DJ,
8      2      IRK, LOC, ITERX, INDEX, I, J, N, N
9      C      CALLS OTHER SUBROUTINES IN A BLOCK
10     C6969    FORMAT(1H 8HSUB RC14)
11      CALL RC3
12      CALL RC8
13      CALL RC9
14      CALL RC4
15      CALL RC5 (CPB, F(4), C(4), BC(2))
16      CALL RC21 (CHB(2), F(3), F(4), C(4), CH(2), CPH(2))
17      CALL RC19 (CPB, CHB(2), CC(2), BC(1), F(4))
18      CALL RC5 (CPT, F(6), C(7), BC(3))
19      CALL RC21 (CHB(3), F(5), F(6), C(7), CH(3), CPH(3))
20      CALL RC19 (CPT, CHB(3), CC(3), BC(1), F(6))
21      CALL RC10
22      CALL RC20
23      CALL RC11
24      RETURN
25      END

```

IPRT,5 RC15

ORIGINAL PAGE IS  
OF POOR QUALITY

JB6-G03432\*TPFS.RC15

```

1      SUBROUTINE RC15
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/2/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DI,
8      2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      C6969      FORMAT(1H 8HSUB RC15)
10     C      NAMELIST/SCH/SV
11     C      SHIFTS VALUES IN SV ARRAY
12     DO 1530 I = 1,18
13     DO 1520 J = 1,49
14     JM = SI - J
15     JMM = JM - 1
16     SV(I,JM) = SV(I,JMM)
17     1520 CONTINUE
18     1530 CONTINUE
19     RETURN
20     END

```

IPRT,S RC16

IB6-G03432-TPFS-RC16

```
1      SUBROUTINE RC16
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/Z/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DT,
8      2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      COMMON/R/ XDS, XMH, CAT, WORK, DUM1, DUM2, DUM3, WORK2, RMTB, RMTB2, TIMEOF
10     1      ,RMLIN
11     C6969      FORMAT(1H 8H SUB RC16)
12     C      SETS VALUES FOR SV ARRAY
13     C      ARTERIAL CO2 CONCENTRATION.
14     SV(1,1) = CC(1)
15     C      ARTERIAL O2 CONCENTRATION.
16     SV(2,1) = F(9)
17     C      BRAIN VENOUS CO2 CONCENTRATION.
18     SV(4,1) = CC(2)
19     C      ARTERIAL N2 CONCENTRATION.
20     SV(3,1) = F(10)
21     C      BRAIN VENOUS O2 CONCENTRATION
22     SV(5,1) = F(12)
23     C      BRAIN VENOUS N2 CONCENTRATION.
24     SV(6,1) = C(6)
25     C      TISSUE VENOUS CO2 CONCENTRATION.
26     SV(7,1) = CC(3)
27     C      TISSUE VENOUS O2 CONCENTRATION.
28     SV(8,1) = F(13)
29     C      TISSUE VENOUS N2 CONCENTRATION.
30     SV(9,1) = C(9)
31     C      CARDIAC OUTPUT.
32     SV(10,1) = C(10)
33     C      CEREBRAL BLOOD FLOW.
34     SV(11,1) = C(11)
35     C      TISSUE BLOOD FLOW.
36     SV(12,1) = QF(1)
37     C      ARTERIAL H+ CONCENTRATION.
38     SV(13,1) = CH(1)
39     C      ARTERIAL O2 TENSION.
40     SV(14,1) = F(1)
41     C      INITIAL TIME.
42     SV(15,1) = 0.0
43     C      TOTAL GAS CONCENTRATIONS AT BRAIN EXIT.
44     SV(16,1) = SV(4,1) + SV(5,1) + SV(6,1)
45     C      TOTAL GAS CONCENTRATIONS AT TISSUE EXIT.
46     SV(17,1) = SV(7,1) + SV(8,1) + SV(9,1)
47     C      SIMULATED TIME.
48     SV(18,1) = C(35)
49     RETURN
50     END
```

IPRT,5 RC17

JB6-G03432\*TPF5.RC17

```

1      SUBROUTINE RC17
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3          1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4          2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5          3      DQ(4)
6      COMMON/2/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7          1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DT,
8          2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      COMMON/R/ XDS, XMH, CAT, WORK, DUM1, DUM2, DUM3, WORK2, RMTB, RMTB2, TIMEOF
10     1      , RMLIN
11     C      NAMELIST/BAD/CH(4), CADK, D(11), C(12), BC(4), C(37), C(38), VTRAN(14),
12     C      1TFRM, VI, C(20), C(16), VTRAN(15), C(21), VTRAN(13), C(37), D(9), C(11),
13     C      2VTRAN(16), QF(1), VTRAN(17), C(10), F(11),
14     C6969      FORMAT(1H 8HSUB RC17)
15     C      CALCULATES VENTILATION
16     C      CFS H+ CONCENTRATION , EQUATION 6.1 .
17     CH(4) = CADK*D(11)*C(12)/BC(4)
18     IF (C(37) .GT. 1.0E-5)      GO TO 1708
19     1704 VI = C(38)
20     GO TO 1730
21     1708 TERM = 0.0
22     C      DECISION ON ARTERIAL O2 TENSION AT CAROTID BODIES*SITE.
23     IF (VTRAN(14) - 104.0)      1710, 1720, 1720
24     1710 TERM = (23.6E-9)*((104.0 - VTRAN(14))**4.9)
25     C      CONTROLLER EQUATION AS A FUNCTION OF HUMORAL TERMS.
26     1720 VI = C(20)*(C(16)*VTRAN(15) + (1.0 - C(16))*CH(4))
27     1      + C(21)*VTRAN(13) + TERM - C(37)
28     C      INCLUSION OF NEURAL COMPONENT AS A FUNCTION OF WORK LOAD.
29     SVNT2=SSVENT(SSO2W(WORK)) -VI
30     IF((SVNT2.GT.0.).AND.(SVNT2.LE.15.)) VI=VI+SVNT2
31     IF(SVNT2.GT.15.) VI=VI+15.
32     C
33     C      DESCRIPTION OF TRANSIENT VENTILATION RESPONSE.
34     SVNT =SSVENT(RMLIN.) -VI
35     IF(SVNT.GT.0.5) VI=VI+0.75*SVNT
36     C
37     C      EXPIRED VENTILATION RATE, EQUATION 11.1 .
38     1730 VE = VI + D(9)*(C(11)*VTRAN(16) + QF(1)*VTRAN(17) - C(10)*F(11))
39     IF (VI .LT. 0.0 .OR. VE .LT. 0.0)      GO TO 1740
40     RETURN
41     1740 VI = 0.0
42     VE = 0.0
43     RETURN
44     END

```

IPRT,S RC19

DB6-G03432\*TPFS,RC19

```

1      SUBROUTINE RC19 (CPA, CVHBA, CVC, BHCA, FC)
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/2/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPR, CPT, CADK, X, DT,
8      2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      C      NAMELIST/DM2/CPA,CVHBA,CVC,BHCA,FC
10     C6969      FORMAT(1H 8HSUB RC19)
11     C      ITERATES FOR VENOUS BRAIN AND VENOUS TISSUE CO2 CONCENTRATION
12     C      TERM USED IN EQUATION 4.2 .
13     1910 X = (CVC - FC)/(0.01*CPA)
14     C      LOGARITHM SUBROUTINE.
15     X = RCF1(X)
16     C      EQUATION 4.2 .
17     X = BHCA + 0.375*(C(17) - CVHBA) - D(8)*(X - 0.14) + FC
18     CALL RC6 (CVC)
19     CVC = CVC + 2.0*(X - CVC)/3.0
20     IF (ITERX)      1920, 1910, 1920
21     1920 CONTINUE
22     RETURN
23     END

```

PRINT, S RC20



ORIGINAL PAGE IS  
OF POOR QUALITY

JB6-G03432\*TPFS-RC20

```

1      SUBROUTINE RC20
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/Z/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DT,
8      2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      C      NAMELIST/NMF/F
10     C6969      FORMAT(1H 8MSUB RC20)
11     C      SETS TIME DEPENDENT EXPRESSIONS
12     C      ARTERIAL OXYGEN CONCENTRATION INCLUDING EFFECTS OF HEMOGLOBIN.
13     F(9) = D(6)*C(2) + CHB(1)
14     C      ARTERIAL NITROGEN CONCENTRATION.
15     F(10) = D(7)*C(3)
16     C      TOTAL ARTERIAL GAS CONCENTRATION AT LUNG EXIT.
17     F(11) = CC(1) + F(9) + F(10)
18     C      VENOUS BRAIN OXYGEN CONCENTRATION INCLUDING EFFECTS OF HEMOGLOBIN.
19     F(12) = C(5) + CHB(2)
20     C      VENOUS TISSUE OXYGEN CONCENTRATION INCLUDING EFFECTS OF HEMOGLOBIN.
21     F(13) = C(8) + CHB(3)
22     C      OXYGEN TENSION IN BRAIN.
23     F(17) = C(5)/D(3)
24     C      NITROGEN TENSION IN BRAIN.
25     F(18) = C(6)/D(4)
26     C      PRODUCT OF DIFFUSION COEFFS. AND DIFFERENTIAL BRAIN - CSF GAS TENSIONS
27     F(14) = C(27)*(CPB - C(12))
28     F(15) = C(28)*(F(17) - C(13))
29     F(16) = C(29)*(F(18) - C(14))
30     C
31     RETURN
32     END

```

SPRT,S RC21

DB6-G03432\*TPFS\*RC21

```

1      SUBROUTINE RC21 (CHBA, FA, FD, CCA, CHA, CPHA)
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/Z/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DT,
8      2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      C6969      FORMAT(1H 8H SUB RC21)
10     C      NAMELIST/PR/CHBA,FA,FD,CCA,CHA,CPHA
11     C      COMPUTES H+ ION, PH, AND OXYHEMOGLOBIN
12     C      ARTERIAL H+ CONCENTRATION.
13     C      CHA = CADK*FD/(CCA - FD)
14     C      ARTERIAL PH.
15     C      CPHA = 9.0 - RCF1(CHA)
16     C      DEVELOPMENT OF EXPRESSION USED IN CALCULATION OF ARTERIAL
17     C      OXYHEMOGLOBIN SATURATION.
18     C      X = RCF2(CPHA)
19     C      X = -X * FA
20     C      X = (1.0 - EXP (X))**2
21     C      X=ABS(X)
22     C
23     C      ARTERIAL OXYHEMOGLOBIN CONCENTRATION.
24     C      CHBA = X*C(17)
25     C      RETURN
26     C      END

```

PRRT,S RCF1

ORIGINAL PAGE IS  
OF POOR QUALITY

DB6-G03432-TPP9-RCF1

```

1      FUNCTION RCF1(W)
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/2/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DT,
8      2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      C      LOGARITHM TO BASE 10
10     RCF1 = 0.43429448 * ALOG(W)
11     RETURN
12     END

```

SPRT,S RCF2

ORIGINAL PAGE IS  
OF POOR QUALITY

DB6-G03432\*TPFS.RCF2

```

1      FUNCTION RCF2(Z)
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/Z/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DT,
8      2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      C      OXYHEMOGLOBIN = PH EMPIRICAL FUNCTION
10     C      EQUATION 3.4 .
11     RCF2 = ((0.0066815*Z) - 0.10098)*Z + 0.44921)*Z - 0.454
12     RETURN
13     END

```

SPRT,S RCF3

ORIGINAL PAGE IS  
OF POOR QUALITY

IB6-GD3432\*TPF5.RCF3

```

1      FUNCTION RCF3(KK)
2      DIMENSION C(40), XN(40,2), SV(18,50), VTRAN(18), RK(14,4),
3      1      SC(14,5), DC(14), A(6), D(15), F(20), VOL(10), RMT(2),
4      2      BC(4), QF(6), TAU(5), CC(3), CHB(3), CH(4), CPH(3),
5      3      DQ(4)
6      COMMON/Z/ C, XN, SV, VTRAN, RK, SC, DC, A, D, F, VOL, RMT, BC, QF,
7      1      TAU, CC, CHB, CH, CPH, DQ, VE, VI, CPB, CPT, CADK, X, DT,
8      2      IRK, LOC, ITERX, INDEX, I, J, M, N
9      C      VTRAN FUNCTION
10     C      VARIABLES WITH TIME DELAYS USED IN EQUATIONS 8.1-8.1 .
11     RCF3 = SV(KK,LOC) + (SV(KK,LOC + 1) - SV(KK,LOC))*DT/D(14)
12     RETURN
13     END

```

IPRT,S SS024

ORIGINAL PAGE IS  
OF POOR QUALITY

DB6-G03432\*TPFS.SS02W

```
1      FUNCTION SS02W(X)
2      C  CALCULATION OF STEADY-STATE OXYGEN REQUIREMENTS FOR VARIOUS LEVELS
3      C  OF WORK LOAD (X=WATTS).
4          IF(X.GT.250.)GOTO1
5          IF( X.LT.75.)GOTO2
6          SS02W=-.072+(X/70.)
7          RETURN
8      1    SS02W=3.5
9          RETURN
10     2    SS02W=(X/75.)+.215*(75.-X)/75.
11     C
12     RETURN
13     END
```

DPRT.S SSVENT

DB6-603432\*TPFS\*SSVENT

```
1      FUNCTION SSVENT(X)
2      C  CALCULATION OF STEADY-STATE VENTILATION RATE AS A FUNCTION OF TISSUE
3      C  OXYGEN METABOLIC RATE.
4      IF(X.LE..215) SSVENT=5.398
5      IF((X.GT..215).AND.(X.LT.2.))SSVENT=25.*X
6      IF(X.GE.2.)SSVENT=50.+50.*(X-2.)
7      C
8      RETURN
9      END
```

DPRT,S PLOT

DB6-G03432\*TPFS.PLOT

```

1      SUBROUTINE PLOT
2      COMMON/PLTBUF/IPLTPT,TLOTBF(181),PLOTBF(181,8),ALPHA(8),KSTOPP
3      & ,TMAX,KPLT
4      DIMENSION HEAD(6)
5      DATA HEAD/'GRODINS RESPIRATORY CONTROL MODEL  '/
6      DIMENSION TSTEPX(4),TSTALP(4)
7      DATA TSTEPX/1.,1.,60.,1440./
8      DATA TSTALP/'MINS','MINS','HOUR','DAYS'/
9      DIMENSION IIDAS(8,2),PPARS(8,3)
10     DIMENSION XNO(5)
11     DATA ISA/'S  '/,IN/'N  '/,IY/'Y  '/,IPLL/'PLOT'/
12     DATA IFIRST,KSTOPP,NOP/0,0,0/
13     IF (IFIRST.EQ.0) CALL INIT7(300)
14     K=8
15     IFIRST=1
16     IEXECI=IPLL
17     IS = IPLTPT - 1
18     IF(IPLTPT.LE.1) GO TO 1
19     IF(IEXECI.EQ.IPLL.AND.KSTOPP.EQ.0)
20     1 GO TO 28
21     IF(IEXECI.EQ.IPLL.AND.KSTOPP.GE.1)
22     1 GO TO 120
23     1 IS = 1
24     11 WRITE(6,2)
25     2 FORMAT(9X,' GRAPHIC OUTPUT(Y,N,S),TIME INTERVALS,STARTX,'
26     & ,',STOPX,(A2,3F5.0),..')
27     4 READ(5,3,ERR=4) I,TSPLOT,STARTX,STOPX
28     3 FORMAT(A2,3F5.0)
29     IF(I.EQ.IN) RETURN
30     IF(I.EQ.ISA) GO TO 10
31     IF (STOPX,LT.1.) STOPX=TMAX
32     IF (STARTX,LT.1.) STARTX=TLOTBF(1)
33     RUNSTP=STOPX-STARTX
34     TTSPT = TSPLOT
35     IF(I.NE.IY) GO TO 11
36     IF(TSPLOT.LT.3..OR.TSPLOT.GT.5.) GO TO 11
37     C   HERE TO BUILD PLOT PARAMETERS.
38     NOP = 0
39     DO 20 I = 1,K
40     17 WRITE(6,5) ALPHA(I)
41     5 FORMAT(' ',A6,13X,'Y SCALE  (A4,8X,F4.0,2F6.0)'/
42     1 ' PLOT(Y,N,S) LOC HIGH LOW  ...')
43     7 READ(5,6,ERR=7) I1,XP1,XP2,XP3
44     6 FORMAT(A4,8X,F4.0,2F6.0)
45     IF(I1.EQ.IY) GO TO 18
46     IF(I1.EQ.ISA) GO TO 19
47     IF(I1.NE.IN) GO TO 17
48     PPARS(I,1) = 0.
49     GO TO 20
50     18 IF(XP1.LE.0.) GO TO 17
51     IF(IFIX(XP1).GT.8) GO TO 17
52     IF(XP3.GE.XP2) GO TO 17
53     PPARS(I,1) = XP1
54     PPARS(I,2) = XP2
55     PPARS(I,3) = XP3
56     19 IF(IFIX(PPARS(I,1)).GT.NOP) NOP = PPARS(I,1)

```



ORIGINAL PAGE IS  
OF POOR QUALITY

```

57      IF(IPARS(I,1) .LT. 1.) GO TO 17
58      20 CONTINUE
59      C   BUILD WHOLE PAGE GRAPH.
60      10 CALL NEWPAG
61          CALL MOVABS(300,775)
62          CALL DMPBUF
63          WRITE (6,21) HEAD
64      21 FORMAT(25X,6A6)
65          CALL MOVABS(3,100)
66          CALL DRWABS(1000,100)
67          CALL DRWABS(1000,750)
68          CALL DRWABS(3,750)
69          CALL DRWABS(3,100)
70          CALL MOVABS(300,750)
71          CALL DRWABS(300,100)
72          CALL DMPBUF
73          I1 = TTSP - 1.
74          I2 = 300
75          I3 = 700 / (I1 + 1)
76          DO 25 I = 1,I1
77              I2 = I2 + I3
78              CALL MOVABS(I2,100)
79              CALL DRWABS(I2,750)
80      25 CONTINUE
81          CALL DMPBUF
82      C   DRAW DIVISIONS FOR DIFFERENT GRAPHS.
83          I3Y = 650 / NOP
84          I4Y = I3Y / 2
85          I1 = NOP - 1
86          I2 = 100
87          IF(I1 .LT. 1) GO TO 280
88          DO 26 I = 1,I1
89              I4 = I2 + I4Y
90              CALL MOVABS(300,I4)
91              CALL DRWABS(305,I4)
92              CALL MOVABS(995,I4)
93              CALL DRWABS(1000,I4)
94              I2 = I2 + I3Y
95              CALL MOVABS(3,I2)
96              CALL DRWABS(1000,I2)
97              CALL DMPBUF
98      26 CONTINUE
99      280 I4 = I2 + I4Y
100          CALL MOVABS(300,I4)
101          CALL DRWABS(305,I4)
102          CALL MOVABS(995,I4)
103          CALL DRWABS(1000,I4)
104          CALL DMPBUF
105      C   ADD ALPHA
106          I1 = 750
107          IDUP = 1
108          CALL VWINDO(0.,1023.,0.,780.)
109          CALL SWINDO(0,1023,0,780)
110          DO 70 I = 1,NOP
111              I11 = I1
112              I22 = I1 - I3Y + 30
113              ID = 0

```

```

114      DO 60 J = 1,K
115      IF(IFIX(PPARS(J,1)) .NE. 1) GO TO 60
116      CALL MOVABS(0,111)
117      CALL ANMODE
118      CALL DMPBUF
119      WRITE(6,55) ALPHA(J),PPARS(J,2)
120      55 FORMAT(' ',6X,A6,F8.2)
121      Y1 = 111 - 11.
122      CALL MOVEA(10.,Y1)
123      IF(ID .EQ. 0) GO TO 410
124      ID1 = 0
125      ID2 = 0
126      CALL DSHARC(104.,Y1,ID,ID1,ID2,IDUP)
127      GO TO 411
128      410 CALL DRAWA(104.,Y1)
129      411 CALL MOVABS(0,122)
130      CALL ANMODE
131      CALL DMPBUF
132      WRITE(6,55) ALPHA(J),PPARS(J,3)
133      Y1 = 122 - 13.
134      CALL MOVEA(10.,Y1)
135      IF(ID .EQ. 0) GO TO 413
136      ID1 = 0
137      ID2 = 0
138      CALL DSHARC(104.,Y1,ID,ID1,ID2,IDUP)
139      GO TO 415
140      413 CALL DRAWA(104.,Y1)
141      415 ID = ID + 1
142      111 = 111 - 21
143      122 = 122 - 21
144      60 CONTINUE
145      11 = 11 - 13Y
146      70 CONTINUE
147      X = (RUNSTP/TTSPT) + .000001
148      I = 1
149      C      IF(X .LT. 1.) I = 1
150      C      IF(X .GE. 60.) I = 3
151      C      IF(X .GE. 1440.) I = 4
152      X = X * TSTEPX(I)
153      IF (STARTX.LT.1.) STARTX=TLOTBF(1)
154      X2 = STARTX / TSTEPX(I)
155      X1 = X2 + X
156      11 = TTSPT
157      DO 75 J = 1,11
158      XNO(J) = X1
159      X1 = X1 + X
160      75 CONTINUE
161      CALL MOVABS(0,100)
162      CALL ANMODE
163      CALL DMPBUF
164      IF(TTSPT .GT. 4.) GO TO 80
165      IF(TTSPT .GT. 3.) GO TO 79
166      WRITE(6,93) X2,(XNO(J),J=1,3)
167      93 FORMAT(' ',15X,F6.1,11X,F6.1,11X,F6.1,10X,F6.1)
168      GO TO 85
169      79 WRITE(6,77) X2,(XNO(J),J=1,4)
170      77 FORMAT(' ',15X,F6.1,6X,F6.1,7X,F6.1,6X,F6.1,7X,F6.1)

```

ORIGINAL PAGE IS  
OF POOR QUALITY

```

171      GO TO 85
172      80 WRITE(6,81) X2,(XN0(J),J=1,5)
173      81 FORMAT(' ',11X,6(4X,F6.1))
174      85 WRITE(6,86) T5TALP(1)
175      86 FORMAT(' ',30X,'TIME ('',A4,'')')
176      C INITIALIZE DASH INFORMATION.
177          DO 90 J = 1,6
178          IIDAS(J,1) = 0
179      90 CONTINUE
180      C      IF(IEXECI .EQ. IPLL) RETURN
181      C PLOT VARIABLES...
182      28 IL = 750
183          IS2 = IS
184          I2 = IPLTPT
185          DO 40 I = 1,NOP
186          IL = IL - I3Y
187          IDAS = -1
188          IDUP = 1
189          DO 30 II = 1,K
190          IF(IFIX(PPARS(II,1)) .NE. 1) GO TO 30
191          IDAS = IDAS + 1
192          X = TLOTBF(1)
193          IF (STARTX.GT.0.1) X=STARTX
194          XL = RUNSTP
195          Y = PPARS(II,3)
196          YL = PPARS(II,2) - Y
197          CALL VWINDO(X,XL,Y,YL)
198          CALL SWINDO(300,700,IL,I3Y)
199          X = TLOTBF(IS)
200          -Y = PLOTBF(IS,II)
201          IF(IEXECI .NE. IPLL) GO TO 97
202          IF(Y .LE. PPARS(II,2)) GO TO 105
203          IF(PLOTBF(IS2,II) .GE. PPARS(II,2)) GO TO 30
204          IF (TLOTBF(II).GT.SIOPX) GO TO 30
205          YY = PPARS(II,2)
206          GO TO 96
207      105 IF(Y .GE. PPARS(II,3)) GO TO 97
208          IF(PLOTBF(IS2,II) .LT. PPARS(II,3)) GO TO 30
209          YY = PPARS(II,3)
210      96 CALL MOVEA(X,YY)
211          IDUP = IDUP + 1
212      97 IDUP = IDUP + 1
213          CALL MOVEA(X,Y)
214          DO 35 III = IS2,I2
215          X = TLOTBF(III)
216          Y = PLOTBF(III,II)
217          IF(IDAS .GT. 0) GO TO 33
218          IDUP = IDUP + 1
219          CALL DRAWA(X,Y)
220          GO TO 34
221      33 L = IIDAS(II,1)
222          N = IIDAS(II,2)
223          CALL DSHARC(X,Y,IDAS,L,N,IDUP)
224          IIDAS(II,1) = L
225          IIDAS(II,2) = N
226      34 IF(IDUP .LT. 13 ) GO TO 35
227          CALL DMPBUF

```

```
228      IDUP = 1
229      35 CONTINUE
230      30 CONTINUE
231      40 CALL DMPBUF
232      IF (IDUP .GT. 1) CALL DMPBUF
233      IF (KSTOPP.NE.1) RETURN
234      120 CALL PAGE3
235      IEXECI = 0
236      IF (KSTOPP.GT.1) IEXECI=1PLL
237      GO TO 1
238      END
```

ORIGINAL PAGE IS  
OF POOR QUALITY

6PRT,S GRODAT

DB6-G03432\*TPFS%GRODAT

1			
2		.17827	
3		.53459	
4		.28714	
5		.64121	
6		.00116	
7		.00105	
8		.61553	
9		.00147	
10		.00105	
11		6.0	
12		.74913	
13		48.17427	
14		36.69498	
15		61.17176	
16		40.0	TMAX
17		0.0	CENT SENS PT
18		0.2	HB
19		0.1	R1
20		0.1	R2
21		1.138	CNT SENS COF
22		1.154	CRD BDY SCF
23		3.0	KL
24		1.0	KB
25		39.0	KT
26		0.05	MRB(CO2)
27		0.05	MRB(O2)
28		81.99	D(CO2)
29		4.361	D(O2)
30		2.524	D(N2)
31		260.0	B
32		.0192	FI(CO2)
33		.70	FI(O2)
34		.2808	FI(N2)
35		.1	KCSF
36		0.	I
37		.0078125	H
38		87.55	VI(N)
39		5.39	VI(SS)
40		.25	PRINT ALL TIM
41		000.	UNKNOWN
42		.547	BHC03 BLOOD
43		.585	BHC03 BRAIN
44		.585	BHC03 TISSUE
45		.585	BHC03 CSF
46		.182	RMT(CO2)
47		.215	RMT(O2)
48		0.0	DJ1
49		0.0	DJ2
50	0.	.5	
51	40.	.25	
52	100.	.25	
53			

ORIGINAL PAGE IS  
OF POOR QUALITY

)B6-G03432\*TPFS.GRODAT

1	TTY		
2		.17827	
3		.53459	
4		.28714	
5		.64121	
6		.00116	
7		.00105	
8		.61553	
9		.00147	
10		.00105	
11		6.0	
12		.74913	
13		48.17427	
14		36.69498	
15		61.17176	
16		40.0	TMAX
17		0.0	CENT SENS P1
18		0.2	HB
19		0.1	R1
20		0.1	R2
21		1.138	CNT SENS C0F
22		1.154	CRTD BDY SCF
23		3.0	KL
24		1.0	KB
25		39.0	K1
26		0.05	MRB(C02)
27		0.05	MRB(02)
28		81.99	D(C02)
29		4.361	D(02)
30		2.524	D(N2)
31		260.0	B
32		.0192	FI(C02)
33		.70	FI(02)
34		.2808	FI(N2)
35		.1	KCSF
36		0.	T
37		.0078125	h
38		87.55	VI(N)
39		5.39	VI(SS)
40		.25	PRINT ALL TIM
41		000.	UNKNOWN
42		.547	BHC03 BLOOD
43		.585	BHC03 BRAIN
44		.585	BHC03 TISSUE
45		.585	BHC03 CSF
46		.182	RMT(C02)
47		.215	RMT(02)
48		0.0	DJ1
49		0.0	DJ2

ORIGINAL PAGE IS  
OF POOR QUALITY

@PRT,S SEADAT

JB6-G0343Z\*TPFS\*SEADAT

1 TTY

ORIGINAL PAGE IS  
OF POOR QUALITY

2 .05269  
 3 .15144  
 4 .79587  
 5 .63977  
 6 .00114  
 7 .00974  
 8 .61323  
 9 .00145  
 10 .00974  
 11 6.00000  
 12 .73723  
 13 47.85777  
 14 36.01092  
 15 567.44715  
 16 40.0  
 17 0.0  
 18 0.2  
 19 0.1  
 20 0.1  
 21 1.138  
 22 1.154  
 23 3.0  
 24 1.0  
 25 39.0  
 26 0.05  
 27 0.05  
 28 81.99  
 29 4.361  
 30 2.524  
 31 760.0  
 32 .0004  
 33 .2096  
 34 .7900  
 35 .1  
 36 0.  
 37 .0078125  
 38 87.55  
 39 5.39  
 40 .25  
 41 000.  
 42 .547  
 43 .585  
 44 .585  
 45 .585  
 46 .182  
 47 .215  
 48 0.0  
 49 0.0

TMAX  
 CENT SENS P1  
 HB  
 R1  
 R2  
 CNT SENS COF  
 CRTD RDY SCF  
 KL  
 KB  
 KT  
 MRB(CO2)  
 MRB(O2)  
 D(CO2)  
 D(O2)  
 D(N2)  
 B  
 FI(CO2)  
 FI(O2)  
 FI(N2)  
 KCSF  
 T  
 H  
 VI(N)  
 VI(SS)  
 PRINT ALL TIM  
 UNKNOWN  
 BHC03 BLOOD  
 BHC03 BRAIN  
 BHC03 TISSUE  
 BHC03 CSF  
 RMT(CO2)  
 RMT(O2)  
 DJ1  
 DJ2

SPRT,S CO2DAT

DB6-G03432\*TPFS.CO2DAT

1	TTY		
2		.05269	
3		.15144	
4		.79587	
5		.63977	
6		.00114	
7		.00974	
8		.61323	
9		.00145	
10		.00974	
11		6.00000	
12		.73723	
13		47.85777	
14		36.01092	
15		567.44715	
16		40.0	TMAX
17		0.0	CENT SENS P1
18		0.2	H8
19		0.1	R1
20		0.1	R2
21		1.138	CNT SENS COF
22		1.154	CRTD BDY SCF
23		3.0	KL
24		1.0	KB
25		39.0	KT
26		0.05	MRB(CO2)
27		0.05	MRB(O2)
28		81.99	D(CO2)
29		4.361	D(O2)
30		2.524	D(N2)
31		760.0	B
32		.0700	FI(CO2)
33		.2096	FI(O2)
34		.7204	FI(N2)
35		.1	KCSF
36		0.	T
37	.0078125		H
38		87.55	VI(N)
39		5.39	VI(SS)
40		.25	PRINT ALL TIM
41		000.	UNKNOWN
42		.547	BHC03 BLOOD
43		.585	BHC03 BRAIN
44		.585	BHC03 TISSUE
45		.585	BHC03 CSF
46		.182	RMT(CO2)
47		.215	RMT(O2)
48		0.0	DJ1
49		0.0	DJ2

ORIGINAL PAGE IS  
OF POOR QUALITY

PRT,S RUNB